www.bego.com

VarseoWax CAD/Cast

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'EMPLOI
INSTRUCCIONES DE USO
ISTRUZIONI PER L'USO

Partners in Progress

VarseoWax CAD/Cast

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
Achtung

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VarseoWax CAD/Cast

VarseoWax CAD/Cast: printing resin for creating burnout 3D-printed partial denture frames.

1. Intended use/Indication

VarseoWax CAD/Cast is a printing resin for creating burnout 3D-printed partial denture frames. It is restricted for use in dental applications and may only be used by dental staff.

RX only

VarseoWax CAD/Cast is one of the system components in the BEGO Varseo 3D print system and has been optimised for use in the Varseo 3D printer. The printing settings can be found in the equipment instructions for use.

2. Contraindications

VarseoWax CAD/Cast is intended exclusively for creating burnout frames within the scope of partial denture construction. VarseoWax CAD/Cast is not approved for any other indications. According to Medical Device Directive 93/42/EEC concerning medical devices, VarseoWax CAD/Cast is **not** considered as a medical device. VarseoWax CAD/Cast must **not** be brought into direct contact with patients (e.g., trial fitting in patient's mouth)!

No contraindications during processing are known. However, if intolerances or allergies to the resin occur, the user should discontinue use of the material. Please contact your attending doctor accordingly.

3. Safety instructions

VarseoWax CAD/Cast is produced and tested according to the most stringent quality standards. In order to ensure optimum further processing, please read the information contained in the instructions for use carefully. Improper use and failure to observe the information can have a detrimental effect on the quality. Nitrile gloves, goggles and a coat must be worn as a means of protection when handling the resin.

4. Side effects and precautions

No side effects resulting from the use of VarseoWax CAD/Cast are known. However, we cannot completely exclude the possibility of personal reactions to individual components in isolated cases. In such cases, the respective user should discontinue use of VarseoWax CAD/Cast.

Information on hazards as per MSDS

- Causes skin irritations
- Can cause allergic skin reactions
- Causes severe eye irritation
- Toxic to aquatic life with long-lasting effects

Safety instructions as per MSDS

- Avoid inhaling aerosol
- Avoid release to the environment
- Wear protective gloves/protective clothing/eye protection/face protection
- In the event of skin irritation/rash: seek medical advice/attention
- In the event of persistent eye irritation: seek medical advice/attention
- Dispose of contents/container in accordance with local and national regulations

5. General information on handling

Delivery

VarseoWax CAD/Cast is supplied in light-tight, sealed bottles. The filling quantities are:

- REF 41000 = 1 kg
- REF 41020 = 0.5 kg

Please check the following points on receipt of the goods:

- Integrity of the bottle/pack
- Quantity
- Shipping documents and designation

Storage

VarseoWax CAD/Cast must be stored in the original sealed bottle at room temperature (approx. 22 °C) in a dark, dry place. It must be ensured that the temperature does not drop below +5 °C and does not exceed +35 °C! The minimum shelf life date printed on the product must be observed. Perfect processing cannot be guaranteed if materials which have exceeded their minimum shelf life date are used.

6. Processing

Please wear protective gloves (nitrile gloves), protective clothing, goggles and/or face protection during processing!

The ideal working temperature range for VarseoWax CAD/Cast is between 20 and 30 °C. The material must be shaken thoroughly for approx. 5 min before being poured into the clean Varseo container. When decanting, make sure that the printing resin is exposed to daylight for as short a period of time as possible.

For further processing – selecting the resin, setting up the print job – as part of the printing process, follow the instructions for use of the Varseo printer. Before starting the printing procedure, the VarseoWax CAD/Cast must be mixed so as to form a homogeneous mixture. Insufficient mixing can lead to deviations in the colour of the printing resin. Varseo enables automatic mixing before each printing procedure, which is started by preselecting the mix mode in the Varseo printer program. Alternatively, mixing can be performed manually by moving the stamp up and down slowly five times.

Note: If the resin is stored in the bowl, always store it with the stamp on the film to prevent the build-up of pigments on the film.

Subsequent processing

On completion of printing, the print objects are released from the build platform by actuating the ejector and using the spatula supplied. The print object should be cleaned in two steps with ethanol using an ultrasonic bath.

1. Clean the print object for 3 min in a reusable ethanol solution in an unheated ultrasonic bath.
2. The precleaned print object must be cleaned thoroughly for 2 min using a fresh ethanol solution. Ideally, this is done by spraying the print object using a spray flask in order to rinse off any remaining resin residue. The entire cleaning process should not take longer than 5 minutes as this could otherwise have a detrimental effect on the print objects. After cleaning, the print object is dried using compressed air, if possible under suction. If there is liquid resin still adhering to the print object, this can be removed without trace by spraying again with ethanol and drying once more.

Finishing

Following this, the support structures are removed. To this end, either a cutting wheel or side cutters can be used. It must be ensured that the printed object is not deformed!

In order to prevent deformation, we recommend post-curing the thoroughly cleaned print objects on a model. This can be done in a standard light polymerisation unit for dental laboratories.

The curing times may vary depending on the lamp output. The final properties of the print object depend on the post-curing process.

Post-curing			
Material	Time [min]	Wave length [nm]	UV-A output [W]
VarseoWax CAD/Cast	10	315-400	72

The values were determined using a device with the following output; four 18W/71 lamps (Dulux L Blue); and four 18W/78 lamps (Dulux blue UV-A). The determined time is based on the UV-A lamp power. The measured post-curing time takes 10 min.

Note: Post-curing without halogen light is recommended as this can result in irregular curing and even colour changes on the print object.

If a corresponding device is not available, the duration of the post-curing can be calculated on base of the overall output of the UV-A lamps.

Intensity/post-curing conversion table				
	Time [min]	Wave length [nm]	UV-A output [W]	Intensity
VarseoWax CAD/Cast	10	315-400	72	43.2 kJ [W x sec = J]

Conversion for equipment with different UV-A lamp outputs

VarseoWax CAD/Cast	20	315-400	36
VarseoWax CAD/Cast	30	315-400	24
VarseoWax CAD/Cast	40	315-400	18

Note: If the print objects are invested immediately after separation of the support structures, post-curing in a light polymerisation unit is not necessary. Post-curing is strongly recommended in the case of extended storage and transport periods!

Storage and transportation of printed objects

Post-curing is strongly recommended prior to extended storage and transport periods! The cured print objects are ideally stored at room temperature away from light and transported in a suitable, light-tight transport box!

Disposal

The cured, separated material (base plate, support structure) can no longer be used. Cured material can be disposed of as domestic waste. Unused resin or ethanol used for cleaning with resin residues must be disposed of via the local waste disposal authority or a hazardous waste collection point stating the safety data sheet.

7. Material properties and scope of delivery

Physical data			
Colour	opaque yellow	Flexural modulus	1,900 MPa
Viscosity	600–900 mPa*s	Shore hardness	81 D
Density at 22 °C	1.08 g/cm ³	Heat resistance	55°C at 1.8 MPa
Flexural strength	60 MPa	Residual ash content	≤ 0.1% at 700 °C

Scope of delivery				
	Contents	Unit	Qty	REF
VarseoWax CAD/Cast	1 kg	bottle	1	41000
VarseoWax CAD/Cast	0.5 kg	bottle	1	41020

8. Equipment

VarseoWax CAD/Cast has been designed for use in the Varseo printer from BEGO Bremer Goldschlägerei Wilhelm-Herbst GmbH & Co KG.

9. Label symbols

Manufacturer

Batch code

Catalogue number

Protect from sunlight


Consult instructions for use

Minimum shelf life

Warning

Storage and transport temperature

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BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG