**REF** 00880/00881/00881P/00882

**CE** 0483 Typ 3, DIN EN ISO 22674 / DIN EN ISO 9693

**DEU Gebrauchsinformation**

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
**REF** 00880/00881/00881P/00882

 0483 Type 3, DIN EN ISO 22674/  
 DIN EN ISO 9693

**ENG Operating Instructions**
 Please read carefully before use  
 Keep for referral purposes

**1. Intended purpose**

NiCrMo - Dental alloy based on nickel for metal/ceramic bonding, type 3 acc. to DIN EN ISO 22674 /DIN EN ISO 9693 for melting and casting.

**2. Indication**

Non-precious alloy for the crown and bridge technique.

**2.1 Contraindication**

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

**3. Fields of application**

FINOBOND SUPERIOR is a nickel-based cast alloy. It is particularly suitable for the production of crowns and bridges, double crowns, implant-supported superstructures, adhesive bridges and for metal ceramics.

**4. Safety precautions**
**4.1 Contraindications and side-effects**

Signs of intolerance to non-precious FINOBOND SUPERIOR alloys are extremely rare if the manufacturer's Instructions for use are adhered to.

If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation to the eyes and/or respiratory tract. Please adhere to warnings and safety data sheets.

**4.2 Interactions**

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient's mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

**4.3 Warnings and precautions**

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly also burs. Processed material may therefore be hot and/or have sharp edges. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Please consult a doctor if you are exposed to individual elements.

 **Note:**

Please observe the relevant safety data sheet!

**5. Processing instructions**
**5.1 Wax-up**

Construction must be in line with dental regulations in order to fulfil mechanical requirements. Use only materials (e.g. wax, plastic) that burn without residue. Detailed information can be found in the processing instructions.

**5.2 Investing**

Investment materials for non-precious alloys with high overall expansion must be used, such as e.g. FINOVEST SPEED C+B (REF 40125).

**5.3 Melting and casting**

Only suitable ceramic-based melting crucibles should be used for FINOBOND SUPERIOR alloys. Care must be taken to ensure that for each alloy separate melting crucibles are used which are not contaminated and/ or damaged. The use of new material is a prerequisite for ensuring clear batch traceability. Suitable melting and casting processes are high- frequency melting with vacuum pressure casting or centrifugal casting machines, arc melting with vacuum pressure casting or open melting. After the breakup of the melt, cast off immediately to avoid overheating of the alloy. Following the casting process, the casting ring is carefully removed and should cool slowly in the air.

**5.4 Work up**

The finishing of the framework and the surfaces should be carried out by means of finely-toothed hard metal burs, ceramic-bonded grinding points or diamond grinding points.

**5.5 Welding**

Suitable welding processes such as laser welding and TIG impulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters of each individual laser machine.

Suitable welding wires are available as filler material for all FINOBOND SUPERIOR alloys e.g. NiCr welding wire 0.50 mm FINOLLOY SUPERIOR (REF 42067).

**5.6 Soldering**

Soldering should be avoided if possible to reduce the number of materials. If, however, it is necessary to solder, please use a suitable solder and flux which are appropriate for the composition and melting temperature of the alloy used., e.g. FINOLLOY SUPERIOR NiCr (REF 42066) in conjunction with FINO DEOXYD SUPERIOR (REF 42088).

**5.7 Porcelain veneering**

FINOBOND SUPERIOR bonding alloys can be veneered with all suitable ceramic masses, e.g. as long as the CTE is suitable for the alloy. The highest recommended firing temperature is 980 °C / 1796 °F. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

**5.8 Reuse**

The alloy may be changed in its composition and properties if reused. Its reuse is therefore not permitted. This is the only way to ensure clear batch traceability.

**6. Storage**

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of the FINOBOND SUPERIOR alloys.

**7. Physical data**

 Elongation limit 340 MPa; tensile strength 580 MPa; Vickers hardness 195 HV10; ultimate elongation 15 %; modulus of elasticity 170 GPa; density 8.2 g/cm<sup>3</sup>; Solidus temperature 1260 °C / 2300 °F; Liquidus temperature 1350 °C / 2462 °F; Casting temperature (recommended) 1410 °C / 2570 °F; coefficient of thermal expansion 14.0 x 10<sup>-6</sup>K<sup>-1</sup> (25-500 °C / 77-932 °F); Weight per casting approx. 6 g.

\* Variations of ± 10 % are possible depending on the batch.

**7.1 Composition in %:**

Ni	Cr	Mo	Si	Fe
61,0	26,0	11,0	1,5	1,2

Additional elements &lt; 1%: Co

This product contains nickel. It should not be used in patients with a known nickel allergy.

Beryllium, lead and cadmium free.

**8. Delivery forms**

 FINOBOND SUPERIOR porcelain-bonding alloy  
 250 g 00880  
 500 g 00881  
 1000 g 00882  
 Trial pack 50 g 00881P

**9. Guarantee**

Thank you for choosing a quality product from Dentaforum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application.

 As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at [www.dl-shop.com](http://www.dl-shop.com), even though you may frequently use the same product.

**9.1 Quality**

FINO ensures faultless quality of its products. These recommendations are based upon Dentaforum's own experiences.

The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, FINO, have no influence on the processing on site.

**9.2 General information**

These Instructions for use contain the main steps for processing the non-precious FINOBOND SUPERIOR alloys.



**10. Disposal**

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.


**11. Symbols explanation**
**11.1 Manufacturing**

	Item numbers
	Lot code
	Shelf life
	non-sterile
	Unit in the packaging
	Manufacturer
	Date of manufacture

**11.2 Application**

	Please read carefully before use Keep for referral purposes
	Only to be used by qualified personnel!
	Do not reuse
	Biohazard

**11.3 Storage**

	Temperature limit
	Keep away from direct sunlight
	Protect against moisture
	Do not use damaged packaging
	top - store upright
	Caution fragile

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