FINOBOND NF SUPERIOR

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00890/00890P/00891/00892

CE 0483

DIN EN ISO 9693/DIN EN ISO 22674

Deutsch

Gebrauchsinformation

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User Information

1. Fields of application

Nickel-free, CoCr-based. Safe processing for all casting protocols. Universally applicable for all dental technological indications. For the crown and bridge technique, secondary parts of combination work, adhesive bridges, implant-supported super-constructions as well as the double-crown technique.

2. Safety precautions

Note:

For processing with rotary instruments use an object extraction.

Wearing of a face and dust mask is recommended.

2.1 Contraindications and side-effects

Do not use in case of a known sensitising towards one of the components of the alloy. If processed according to user information intolerance occurrences towards this alloy are extremely rare.

In case different groups of alloys are used galvanic effects may occur. Electro-chemically induced irritations have been described for individual cases.

3. Processing instructions

3.1 Wax-up

Varnish the dies with a die spacer before wax-up.

To guarantee a safe spread of the wax-up do not go below a crown wall thickness of 0.5 mm.

3.2 Placement of sprues

For restorations with more than 4 units place sprues in the indirect method.

Recommended length/thickness of wax wires:

Connection between casting canal and core: \varnothing wax wire 3.5 to 4.0 mm.

Cross bar: \varnothing wax wire 4 to 5 mm.

Connector between cross bar and object: \varnothing wax wire 2.5 to 3 mm, length 2.5 to 3 mm. Single crowns and small bridges can be connected directly.

Single crowns: \varnothing wax wire 3 mm, length 15 to 20 mm.

Bigger crowns and bridge pontics: \varnothing wax wire 3.5 to 4.0 mm, length 15 to 20 mm.

3.3 Investing

Any phosphate-bonded graphite-free investment material suitable in the expansion range for non-precious crown and bridge alloys, e.g. FINOVEST SPEED C+B investment material (item-no. 40125), can be used. In case metal mould rings are used a sufficient mould lining is required.

Follow the user information of the investment material.

3.4 Pre-heating

Burning-out of the wax and heating of the mould is done according to the user information of the investment material.

Pre-heating time depends on the mould size and the number of moulds in the pre-heating furnace. Pre-heating temperature 950 °C, holding at final temperature between 30 and 60 min.

3.5 Quantity of metal

Calculation of the required metal quantity according to the following rule of thumb:

Wax weight of the wax-up x density of the alloy (see. 5. Physical data) + 10 g.

3.6 Melting and casting

The alloy can be melted in the high-frequency melting process with vacuum pressure or centrifugal casting units, with an open flame as well as electric arc melting with vacuum pressure casting.

Use only ceramic crucibles and pre-heat them in the furnace (without metal). Use crucibles only for one alloy and clean after each cast (remove casting cinder).

High frequency melting with vacuum pressure casting units:

Adjust vacuum to 250 to 450 mbar. Recommended pressure 450 mbar.

Melt casting cylinders to red heat (no visible beveling of the edges yet). Heat main melt until melt tears and cast immediately. Observe melting process and initiate casting manually.

High frequency melting with centrifugal casting units:

After the last casting cylinder has collapsed wait until the oxide skin tears and cast immediately.

Open flame melting:

Use a shower head on the torch and adjust in a way that the heat is evenly distributed.

If melting is done with propane/oxygen, select maximum output according to the user information of the manufacturer.

If melting is done with acetylene/oxygen, follow the user information of the manufacturer. Acetylene: 0.7 bar/oxygen 3 bar. The length of the blue core of the flame directly at the shower head should be 4 to 5 mm.

During melting after the casting cylinders have collapsed an oxide skin is formed. Melt the metal further until it is visibly moved by the flame pressure under the oxide skin. Start casting process before oxide skin tears.

Avoid overheating of the melt for all casting and melting protocols. Overheating might result in blowholes, coarse grit and micro-porositites.

3.7 Processing

For surface processing use TC cutters. Employ abrasives only for one alloy. Pay attention to a uniform direction of cutting and a low working pressure. Avoid overlapping. Observe maximum allowed speed. Do not go below a minimum cap thickness of 0.2 to 0.3 mm.

3.8 Oxide firing

As an option an oxide firing for the control and cleaning of the cast surface can be done. If not otherwise indicated by the manufacturer of the porcelain this is done at the respective opaque firing temperature with a holding time of 5 min without vacuum. Subsequently blast oxide off with Al₂O₃ (grain approx. 110 µm, e.g. FINOX corundum, item-no. 41033) at a pressure of 2 to 3 bar.

3.9 Soldering

For soldering use FINOLLOY SUPERIOR CoCr solder (item-no. 42065) as well as FINO DEOXYD SUPERIOR flux (item-no. 42068).

3.10 Welding

The alloy is suitable for laser-welding. Suitable wire: FINOLLOY SUPERIOR laser wire (item-no. 42067).

The alloy is also suitable for WIG impulse-welding.

3.11 Porcelain veneering

After processing blast extensively with Al₂O₃ (grain approx. 110 µm, e.g. FINOX corundum, item-no. 41033) at a blasting pressure of 2 to 3 bar. The rough surface improves the bond between alloy and porcelain. The firing of the porcelain materials is done according to the user information of the respective manufacturer. Generally, all suitable porcelains adjusted to the CTE of the alloy can be used for veneering.

3.12 Cooling after porcelain firing

Observe the recommendations of the manufacturer of the porcelain.

3.13 Reuse of old material

Do not use old material when casting the alloy.

Cinder or inclusion of foreign objects (e.g. investment material) may cause defects like porosities or blowholes in the cast structure. Defects result in a reduced hardness of the alloy.

4. Storage

No special measures required.

5. Physical data

Elongation limit 760 N/mm²; tensile strength 940 N/mm²; Vickers hardness 340 HV10; ultimate elongation 5.3 %; modulus of elasticity 211,000 N/mm²; density 8.6 g/cm³; melting interval 1329-1385 °C; CTE 14.4 x 10⁻⁶K⁻¹ (25-500 °C).

5.1 Composition in %:

Co 58.0; Cr 30.0; W 5.5; Mo 3.0; Si 1.5; Ta 1.25.

6. Delivery forms

FINOBOND NF SUPERIOR
porcelain-bonding alloy

	250 g	00890
	500 g	00891
	1000 g	00892
Trial pack	50 g	00890P

7. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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