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**User Information**

**1. Fields of application**

Nickel- and beryllium-free. Specifically developed for veneering with low-fusing porcelains. Applicable in the crown and bridge technique, combination prosthetics as well as the double crown technique.

**2. Safety precautions**

 **Note:**  
For processing with rotary instruments use an object extraction.  
Wearing of a face and dust mask is recommended.

**2.1 Contraindications and side-effects**  
Do not use in case of a known sensitising towards one of the components of the alloy. If processed according to user information intolerance occurrences towards this alloy are extremely rare.  
In case different groups of alloys are used galvanic effects may occur. Electro-chemically induced irritations have been described for individual cases.

**3. Processing instructions**

**3.1 Wax-up**  
Varnish the dies with a die spacer before wax-up.  
To guarantee a safe spread of the wax-up do not go below a crown wall thickness of 0.5 mm.

**3.2 Placement of sprues**  
For restorations with more than 4 units place sprues in the indirect method.  
Recommended length/thickness of wax wires:  
Connection between casting canal and core:  $\varnothing$  wax wire 3.5 to 4.0 mm.  
Cross bar:  $\varnothing$  wax wire 4 to 5 mm.  
Connector between cross bar and object:  $\varnothing$  wax wire 2.5 to 3 mm, length 2.5 to 3 mm.  
Single crowns and small bridges can be connected directly.

Single crowns:  $\varnothing$  wax wire 3 mm, length 15 to 20 mm.  
Bigger crowns and bridge pontics:  $\varnothing$  wax wire 3.5 to 4.0 mm, length 15 to 20 mm.

**3.3 Investing**  
Any phosphate-bonded graphite-free investment material suitable in the expansion range for non-precious crown and bridge alloys, e.g. FINOVEST SPEED C+B investment material (item-no. 40125), can be used. In case metal mould rings are used a sufficient mould lining is required.  
Follow the user information of the investment material.

**3.4 Pre-heating**  
Burning-out of the wax and heating of the mould is done according to the user information of the investment material.  
Pre-heating time depends on the mould size and the number of moulds in the pre-heating furnace. Pre-heating temperature 950 °C, holding at final temperature between 30 and 60 min.

**3.5 Quantity of metal**  
Calculation of the required metal quantity according to the following rule of thumb:  
Wax weight of the wax-up x density of the alloy (see. 5. Physical data) + 10 g.

**3.6 Melting and casting**  
The alloy can be melted in the high-frequency melting process with vacuum pressure or centrifugal casting units, with an open flame as well as electric arc melting with vacuum pressure casting.  
Use only ceramic crucibles and pre-heat them in the furnace (without metal). Use crucibles only for one alloy and clean after each cast (remove casting cinder).  
High frequency melting with vacuum pressure casting units:  
Adjust vacuum to 250 to 450 mbar. Recommended pressure 450 mbar. After glow shade disappears melt for another 2-4 s without tearing the melt. Now initiate the casting process.  
In case of high-frequency melting in centrifugal casting units the casting process is initiated after the glow shade disappears.  
Open flame melting:  
Use a shower head on the torch and adjust in a way that the heat is evenly distributed.  
If melting is done with propane/oxygen, select maximum output according to the user information of the manufacturer. Oxygen: 2 to 3 bar.  
If melting is done with acetylene/oxygen, follow the user information of the manufacturer. Acetylene: 0.7 bar/oxygen 3 bar. The length of the blue core of the flame directly at the shower head should be 4 to 5 mm.

During melting after the casting cylinders have collapsed an oxide skin is formed. Melt the metal further until it is visibly moved by the flame pressure under the oxide skin. Start casting process before oxide skin tears.  
Avoid overheating of the melt for all casting and melting protocols. Overheating might result in blowholes, coarse grit and micro-porosities.

**3.7 Processing**  
For surface processing use TC cutters. Employ abrasives only for one alloy. Pay attention to a uniform direction of cutting and a low working pressure. Avoid overlapping. Observe maximum allowed speed.  
Do not go below a minimum cap thickness of 0.2 to 0.3 mm.

**3.8 Oxide firing**  
As an option an oxide firing for the control and cleaning of the cast surface can be done. If not otherwise indicated by the manufacturer of the porcelain this is done at the respective opaque firing temperature with a holding time of 5 min without vacuum. Subsequently blast oxide off with Al<sub>2</sub>O<sub>3</sub> (grain approx. 110 µm, e.g. FINOX corundum, item-no. 41033) at a pressure of 2 to 3 bar.

**3.9 Soldering**  
For soldering use FINOLLOY SUPERIOR CoCr solder (item-no. 42065) as well as FINO DEOXYD SUPERIOR flux (item-no. 42068).  
Let objects that have been soldered after firing of the porcelain cool slowly. Soldering after porcelain firing is not recommended.

**3.10 Welding**  
The alloy is suitable for laser-welding. Suitable wire: FINOLLOY laser wire (item-no. 42005).  
The alloy is also suitable for WIG impulse-welding.

**3.11 Porcelain veneering**  
After processing blast extensively with Al<sub>2</sub>O<sub>3</sub> (grain approx. 110 µm, e.g. FINOX corundum, item-no. 41033) at a basting pressure of 2 to 3 bar. The rough surface improves the bond between alloy and porcelain. The firing of the porcelain materials is done according to the user information of the respective manufacturer. Generally all suitable porcelains adjusted to the CTE of the alloy can be used for veneering.

**3.12 Cooling after porcelain firing**  
Observe the recommendations of the manufacturer of the porcelain.

**3.13 Reuse of old material**  
Do not use old material when casting the alloy.  
Cinder or inclusion of foreign objects (e.g. investment material) may cause defects like porosities or blowholes in the cast structure. Defects result in a reduced hardness of the alloy.

**4. Storage**  
No special measures required.

**5. Physical data**  
Elongation limit 660 N/mm<sup>2</sup>; tensile strength 950 N/mm<sup>2</sup>; Vickers hardness 315 HV10; ultimate elongation 11 %; modulus of elasticity 200,000 N/mm<sup>2</sup>; density 8.2 g/cm<sup>3</sup>; melting interval 1280-1350 °C; coefficient of thermal expansion 15.9 x 10<sup>-6</sup>K<sup>-1</sup> (25-500 °C).

**5.1 Composition in %:**  
Co 33.0; Cr 30.0; Fe 29.0; Mo 5.0; Mn 1.5.

**6. Delivery forms**  
FINOBOND LF SUPERIOR  
porcelain-bonding alloy

	250 g	00900
	1000 g	00902
Trial pack	50 g	00900P

**7. Guarantee**  
Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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