FINOLLOY SUPERIOR

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00940/00940P/00942

CE 0483

DIN EN ISO 22674 TYP 5

Deutsch

Gebrauchsinformation

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English

User Information

Please read carefully before use
Keep for referral purposes

1. Fields of application

CoCrMoTa-based, carbon-free, extra spring-hard. Can excellently be laser-welded and can be applied in the whole range of partial denture and combined technique. Particularly suitable for rigid constructions requiring high strength.

2. Safety precautions

Note:

For processing with rotary instruments use an object extraction.

Wearing of a face and dust mask is recommended.

2.2 Contraindications and side-effects

Do not use in case of a known sensitising towards one of the components of the alloy. If processed according to user information intolerance occurrences towards this alloy are extremely rare.

In case different groups of alloys are used galvanic effects may occur. Electro-chemically induced irritations have been described for individual cases.

3. Processing instructions

3.1 Wax-up

A frame thickness of 0.4 mm is sufficient for full plates.

The thickness of the wax-up of the base depends on the delicacy of the major connectors. The more delicate they are the thicker the base has to be. The recommendation for the cross-section of sublingual bars is 2 x 4 mm.

Use round wires with a \varnothing of 3.5 to 4.0 mm as sprues. Minimum length 20 mm. In case of very bulky parts it is recommended to apply sprues of a \varnothing of 4.0 mm. Use 2 to 4 sprues for upper restorations and 2 to 3 sprues for lower ones. To provide guidance for the desired direction of flow sprues must be connected to the wax-up in a bent shape. It is possible to cast through the model or from above.

3.2 Investing

Any phosphate-bonded investment material for partial denture alloys, e.g. FINOVEST investment material (item-no. 40136), can be used.

Follow the user information of the investment material.

3.3 Pre-heating

Burning-out of the wax and heating of the mould is done according to the user information of the investment material.

Depending on the investment material, the casting unit and the restoration the pre-heating temperature is between 850 and 1050 °C.

3.4 Quantity of metal

For one partial denture approx. 4 to 5 casting cylinders of 6 g metal each are needed.

3.5 Melting and casting

The alloy can be melted in the high-frequency melting process with vacuum pressure or centrifugal casting units, with an open flame as well as electric arc melting with vacuum pressure casting.

Use only ceramic crucibles and pre-heat them in the furnace (without metal). Use crucibles only for one alloy and clean after each cast (remove casting cinder).

High frequency melting with vacuum pressure casting units:

Adjust vacuum to 250 to 450 mbar. Recommended pressure 450 mbar.

Melt casting cylinders to red heat (no visible bevelling of the edges yet). After tearing of the oxide skin (melt shine) carry on melting for 1-2 s, then initiate the casting process.

In case of high-frequency melting with centrifugal casting units the casting process is initiated 1-2 s after the tearing of the oxide skin (melt shine). Pre-melting of the metal is recommended if larger quantities of casting metal is used.

Open flame melting:

Use a shower head on the torch and adjust in a way that the heat is evenly distributed.

If melting is done with propane/oxygen, select maximum output according to the user information of the manufacturer.

If melting is done with acetylene/oxygen, follow the user information of the manufacturer. Acetylene: 0.7 bar/oxygen 3 bar. The length of the blue core of the flame directly at the shower head should be 4 to 5 mm.

During melting after the casting cylinders have collapsed an oxide skin is formed. Melt the metal further until it is visibly moved by the flame pressure under the oxide skin. Start casting process before oxide skin tears.

Avoid overheating of the melt for all casting and melting protocols. Overheating might result in blowholes, coarse grit and micro-prosities.

3.6 Blasting

Blast critical areas like inner sides of clasps or retentions particularly gently with Al₂O₃ (grain approx. 125 to 250 μ m, e.g. FINOX corundum).

3.7 Processing

Use cross-cut TC cutters for initial rough processing. Use sinter diamonds and ceramic-bonded abrasives for intermediate and final processing. Shining, polishing and finishing. Subsequently clean thoroughly.

3.8 Soldering

For soldering use FINOLLOY CoCr solder (item-no. 42009) as well as FINO DEOXYD flux (item-no. 42060).

3.9 Welding

The alloy is suitable for laser-welding. Suitable wire: FINOLLOY laser wire (item-no. 42005).

The alloy is also suitable for WIG impulse-welding.

3.10 Cleaning

All common denture cleaning agents for partial dentures made from this alloy can be used. Follow the user information of the denture cleaning agent.

3.11 Reuse of old material

For casting of the alloy only use old material once. Use a minimum of 2/3 new material.

4. Storage

No special measures required.

5. Physical data

Elongation limit 700 N/mm²; tensile strength 910 N/mm²; Vickers hardness 340 HV10; ultimate elongation 5 %; modulus of elasticity 210,000 N/mm²; density 8.2 g/cm³; melting interval 1360-1410 °C.

5.1 Composition in %:

Co 62.0; Cr 30.0; Mo 5.5; Si 1.2; Ta 1.2; N < 1.0.

6. Delivery forms

| | |
|---|--------------|
| FINOLLOY SUPERIOR partial denture alloy | |
| | 250 g 00940 |
| | 1000 g 00942 |
| Trial pack | 50 g 00940P |

7. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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