Wironit®

Wironit® extrahart

	Wironit® Seite/page	Wironit® extrahart Seite/page
de Gebrauchsanweisung	2	30
en Instruction for use	3	31
fr Notice d'utilisation	4	32
es Instrucciones de uso	5	33
it Istruzioni per l'uso	6	34
nl Gebruiksaanwijzing	7	35
pt Instruções de utilização	8	36
ru Инструкция по применению	9	37
ar دليل الاستخدام	10	38
bg Инструкция за употреба	11	39
cs Návod k použití	12	40
da Brugsanvisning	13	41
el Οδηγίες χρήσης	14	42
et Kasutusjuhend	15	43
hr Upute za uporabu	16	44
hu Használati utasítás	17	45
ja 使用説明書	18	46
ko 사용 설명서	19	47
lt Naudojimo instrukcija	20	48
lv Lietošanas pamācība	21	49
pl Instrukcja używania	22	50
ro Instrucțiuni de utilizare	23	51
sk Návod na použitie	24	52
sl Navodilo za uporabo	25	53
sv Bruksanvisning	26	54
th คำแนะนำการใช้งาน	27	55
tr Kullanma talimatı	28	56
zh 使用说明	29	57

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Instruction for use

Dental Co-based casting alloy, Type 5 Wironit® is available as cylinders. Wironit® complies with ISO 22674. REF 50030 – 1000 g; REF 50020 – 250 g; REF 50019 – 1100 g; REF 50032 – 24 g sample	
Alloy characteristics	
According to ISO 22674 free of nickel, cadmium, beryllium and lead	
Type (accord. to ISO 22674)	5
Preheating temperature	°C 950–1050
Solidus, liquidus temperature	°C 1260, 1400
Density	g/cm ³ 8.3
Casting temperature	°C 1460
Young's modulus	GPa 220
Proof strength (R _{p0.2})	MPa 520
Elongation after fracture (A ₅)	% 7
Vickers hardness	(HV10) 360
BEGO color code (white)	8
Investment material: phosphate bonded, e. g. Wirovest (REF 51046)	
Crucible material: ceramic	
Veneering ceramic: not veneerable with ceramic	
Flux: e. g. Minoxid (REF 52530)	
Brazing material: Kobalt-Chrom-Lot (REF 52520)	
Laser wire: Wiroweld (REF 50003, 50005)	
Melting powder: Wiro melt (REF 52526)	

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Soldering / brazing: Fix the parts with soldering investment material (e. g. Bellatherm®). The prepared gap shall not exceed 0.2 mm with parallel walls. Use a suitable BEGO flux material. The flux residues and oxides must be etched off. Clean surface thoroughly by steam cleaning or boiling in distilled water.

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Instruction for use

Dental Co-based casting alloy, Type 5 Wironit® extrahart is available as cylinders. Wironit® extrahart complies with ISO 22674. REF 50060 – 1000 g; REF 50050 – 250 g; REF 50045 – 1100 g; REF 50042 – 24 g sample	
Alloy characteristics	
According to ISO 22674 free of nickel, cadmium, beryllium and lead	
Type (accord. to ISO 22674)	5
Preheating temperature	°C 950–1050
Solidus, liquidus temperature	°C 1270, 1400
Density	g/cm ³ 8.2
Casting temperature	°C 1420
Young's modulus	GPa 238
Proof strength (R _{p0.2})	MPa 515
Elongation after fracture (A ₅)	% 7
Vickers hardness	(HV10) 385
BEGO color code (white)	8
Investment material: phosphate bonded, e. g. Wirovest (REF 51046)	
Crucible material: ceramic	
Veneering ceramic: not veneerable with ceramic	
Flux: e. g. Minoxid (REF 52530)	
Brazing material: Kobalt-Chrom-Lot (REF 52520)	
Laser wire: Wiroweld (REF 50003, 50005)	
Melting powder: Wiromelt (REF 52526)	

Intended Use: Partial denture alloys are intended for casting of dental restorations such as partial dentures and combination works.

Indication: Partial denture alloys are intended to treat the condition of missing hard tissue (teeth).

Contraindications: No contraindications are known. However, unwanted biological reactions such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur. In case of known incompatibilities and allergies to contents of the metallic material it should not be used.

Clinical benefit: Artificial replacement of hard tissue (teeth), to restore masticatory functionality (aesthetic and functional).

Warnings: Metal dust is harmful to your health. Grinding and blasting should be done under a suitable vacuum. Respiratory protection of type FFP3-EN149 is recommended!

Precautions: Electrochemically induced sensitivity may very rarely occur when the restoration makes occlusal or proximal contact with a different alloy. This product has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. Safety in the MR environment is unknown. Scanning a patient with this device may result in patient injury.

Patient group: Objects made of the alloy can be used regardless of patient age. The alloy should not be used in cases of known incompatibilities or known allergies to alloy components.

Adverse reactions: No adverse reactions are known. Nevertheless, rare cases of individual reactions to single components can never be excluded completely. In this case, the application should not be continued.

Prescription device: Caution: US Federal law restricts this device to sale by or on the order of a licensed dentist.

Wax up: Always place sprues in the most solid wax-up areas, e. g. at the transition between saddle and base. Provide solid places which the melt can only reach through a thinly modelled area with an additional sprue (Ø 3 mm). In case of bruxism stronger modellation is required. Do not taper the spruing.

Investing: Use only phosphate bonded partial denture investment materials.

Melting / casting: Do not overheat alloy. Use only clean ceramic crucibles and one crucible per alloy. To enable an exact batch traceability, cast new metal only. If applicable use melting powder. Follow the casting device manufacturer's instructions for parameters and casting procedures. After casting, cool the mould slowly.

Grinding: Blast with Korox® 250 at approx. 4 bar. Critical areas, e.g. inner clasp sides and stress breakers, are to be blasted extremely carefully (e.g. Korox® 50 blasting material). Use fine carbide, ceramically bonded stones or sintered diamond milling tools for finishing.

Polishing: Polishing (Eltropol polishing unit, Wirolyt polishing liquid), rubber-polishing (BEGO rubber polisher, black) and finish-polishing (BEGO cobalt chrome polishing paste, blue). Clean surface thoroughly by steam cleaning or boiling in distilled water.

Ceramic veneering: Partial dentures alloys are not intended for veneering with ceramic.

Acrylic veneering: For veneering with acrylic material follow the recommendations of the manufacturer. Clean surface thoroughly by steam cleaning or boiling in distilled water.

Soldering / brazing: Fix the parts with soldering investment material (e. g. Bellatherm®). The prepared gap shall not exceed 0.2 mm with parallel walls. Use a suitable BEGO flux material. The flux residues and oxides must etched off. Clean surface thoroughly by steam cleaning or boiling in distilled water.

Laser welding: If possible use X-seam and filler material. Follow manufacturer's instructions for use and hazard notes of the laser welding devices! Clean surface thoroughly by steam cleaning or boiling in distilled water.

Cleaning/Disinfection: Finalized casted partial dentures should be cleaned using an ultrasonic bath, by steam cleaning or boiling in distilled water.

Limit of Liability: Except where prohibited by law, BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Storage conditions: No special storage conditions are known.

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can be considered as standard values. Our products are subject to a constant further development; therefore, alterations in construction and composition are reserved. Any serious incident that has occurred in relation to the products should be reported to BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG and the competent authority.

US Labeling requirements: The device labeling meets the recommendations of FDA applicable guidance documents.

Instructions for disposal: Disposal procedures as follows:

Device: The assignment of a waste key identification number as per the European Waste Catalogue Ordinance (AVV) must be carried out in consultation with the regional waste disposal contractor. Do not dispose of with household waste.

Packaging: Packaging must be fully emptied and properly disposed of in compliance with statutory regulations. Packaging that is not fully emptied must be disposed in coordination with the regional waste disposal contractor.

Consult instructions for use

Caution

Use-by-date

Batch code

Contains hazardous substance

Non-sterile

Rx only
Only for technical
personnel!

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