cre-active®/cre-active® Gingiva

CE 0197

cre-active®/cre-active® Gingiva

D Gebrauchsinformation

GB User Instructions

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E Modo de empleo

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S Bruksanvisning

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Heraeus Kulzer GmbH
Grüner Weg 11
63450 Hanau (Germany)

info.lab@heraeus.com
www.heraeus-kulzer.com

GB User instructions cre-active®/cre-active® Gingiva

Photocuring cre-active mass for individual characterization of c&b veneers.

Composition of cre-active
Multifunctional methacrylic acid ester, highly disperse silicon dioxide, cadmium-free pigments

Composition of transpafluids
Dimethacrylate, silicon dioxide

Indications
Individual characterizations of c&b veneers, partial and full dentures and implantation work

– enamel cracks, calcification spots, filling imitations, enamel wedges, colored intermediate layer, abraded surfaces, fissures, gingival segments

cre-active
white, polar, umbra, corn, mango, caramelpink, red, maroon, black

cre-active Gingiva
pink, shade 200, R50, light pink, clear

transpafluids
TF1 – high viscosity, clear
TF2 – low viscosity, clear

Processing
The ready-to-use cre-active mass can be applied directly onto the polymerized veneer material. Instead of using the margin component for thin veneer sections, cre-active can be individually mixed and applied directly to the opaque layer. The transpafluids can be mixed in to reduce the color intensity of cre-active. Use pure transpafluid if highly transparent effects are desired. Always apply a layer of veneering material after application of cre-active, T1 or T2. To repair a destroyed dispersion layer, roughen the surface with a diamond burr and brush on a thin but very moist layer of c&b liquid so that it extends beyond the site of repair. Allow 30 s for c&b liquid to react, do not polymerize separately.

Device	Max. curable layer thickness	Polymerization time
Heraflash	0,2 mm	90 s
Dentacolor® XS	0,2 mm	90 s
UniXS®	0,2 mm	90 s

Note:
cre-active must be applied in thin layers to cured veneer material. cre-active mass may only be mixed with cre-active mass/transpafluids. Cured cre-active layers must be covered with veneering material.

Storage
Do not store material above 25°C. Avoid exposure to direct sunlight. Close all containers after use. Do not use the material after the expiry date. Store out of reach of children.

Risks, safety precautions
Irritates the eyes, respiratory system and skin. Keep containers closed tightly and store in a cool, well-ventilated place. In case of contact with the eyes, rinse immediately with plenty of water and consult a physician. In case of contact with skin, wash immediately with plenty of soap and water.

See Signum® processing instructions for further information.

cre-active
Matrix: Multifunctional methacrylic acid ester
Filler substances: 6% (w/w), silicon dioxide nano-range, silanized, photoinitiators, stabilizers, inorganic pigments, some of which with fluorescent properties (cadmium-free)

cre-active Gingiva
Matrix: Multifunctional methacrylic acid ester
Filler substances: 8% (w/w), silicon dioxide nano-range, silanized, bead polymer, photoinitiators, stabilizers, inorganic pigments (cadmium-free)

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cre-active

White, polar, umbra, corn, mango, caramelpink, red, maroon, black.

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cre-active

white, polar, umbra, corn, mango, caramelpink, red, maroon, black.

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cre-active

white, polar, umbra, corn, mango, caramelpink, red, maroon, black

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® = Registered trademark of Heraeus Kulzer GmbH