REF 02149

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REF 02149

ENG **Instruction for use**

Please read the entire document before using the Plastics.

1. Intended purpose

Light-curing, ready-to-use modelling resin for the mould casting technology. Based on polymethyl methacrylate.

2. Indication

For modelling inlays, onlays, crowns, bridges, pin abutments and abutments, as well as for luting soldered objects.

2.1 Contraindication

Do not use if incompatibilities or allergies to constituents have been confirmed.

3. Fields of application

Only for processing in the dental technology according to the intended use.

4. Safety precautions

No particular hazards worth mentioning.

5. Processing instructions

5.1 Preparation

Isolate plaster surfaces with alginate insulation, isolate polyurethane model materials with petroleum jelly. No isolation is required on polished metal or ceramic surfaces.

5.2 Application

Remove modelling resin directly from the rotary piston syringe and apply with a modelling instrument. Curing is performed with commercially available light curing devices. After polymerisation, the modelled object can be processed immediately. The layer thickness per polymerisation must not exceed more than 3 mm.

The procedure for embedding and preheating the investment ring is analogue to the wax technique. The modelling resin burns at a temperature of 480°C without residue. For best casting results, completely remove the oxygen inhibition layer with isopropanol prior to embedding. For soldered objects, connect the solder joints with modelling resin and then fabricate the solder model (solder block). Expel the resin with an open flame.

6. Storage

P410 Protect from sunlight. P411+P235 Store at temperatures not exceeding 5-25 °C. Keep cool.

⚠ Close container tightly immediately after use.

7. Physical data

Colour: blue-green;
Light polymerisation: 320-500 nm;
Polymerisation time: UV/UV-A: approx. 5 min.,
Stroboscopic lamps: approx. 1 min. 30 s;
LED: approx. 20-40 s.

Burns without residue

8. Delivery forms

Packaging units for Plastics are:

- Modelling Resin
- Instruction for use

FINO RESIN LC Modelling Resin 3 g 02149

Modelling Resin are shown in the illustration:

Plastics are supplied without accessories.

9. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development.

We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

10. Exclusion of liability

All rights reserved. This document may not be reproduced, transmitted, transferred, distributed, modified, summarised or translated into any other language graphically, electronically or mechanically, including but not limited to photocopying, recording, information storage and processing in systems, without the prior written consent of FINO GmbH. Copies of the system software are not permitted.

11. Disposal

P501 Dispose of contents/container in accordance with local/regional/national/ international regulations.

11.1 EWC waste code product

200139 MUNICIPAL WASTES (HOUSEHOLD WASTE AND SIMILAR COMMERCIAL, INDUSTRIAL AND INSTITUTIONAL WASTES) INCLUDING SEPARATELY COLLECTED FRACTIONS; separately collected fractions (except 15 01); Plastics.

12. Symbols explanation

12.1 Manufacturing

REF Item numbers
LOT Lot code
 Shelf life
 non-sterile
 Unit in the packaging
 Manufacturer


12.2 Application

 General danger symbol warns of the risk of serious injury in conjunction with the warning words CAUTION, WARNING or DANGER

 Notes
 Important information

 Please read carefully before use
 Keep for referral purposes
 Do not reuse

12.3 Storage

 Temperature limit
 Keep away from direct sunlight
 Protect against moisture
 top - store upright
 Caution fragile

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