
Attachment bond®


CE 0197

Attachment bond®

- D** Gebrauchsinformation
- GB** User Instructions
- F** Mode d'emploi
- E** Modo de empleo
- I** Istruzioni di impiego
- P** Instruções de uso
- NL** Gebruiksaanwijzing
- S** Bruksanvisning
- DK** Brugsanvisning
- GR** Οδηγίες χρήσεως
- RUS** Информация по использованию

Passung Stahlplatte/Fit to cast framework/Ajustement de la plaque d'acier/Precision della piastra metallica/Ajuste plaqueta de acero/Ajuste da estrutura metálica/Passing stalen plaat/Passing metallskelett/Metallsatulan sovitus/Εφαρμογή σε χυτό οκέτο/Принасовка стальной пластины

Passung Sekundärteil/Fit of attachment/Ajustement de la pièce secondaire/Area di giunzione dell'elemento secondario/Ajuste partes secundarias/Ajuste do Attachment/Passing secondary onderdeel/Passing attachment/Sekundariosan sovitus/Εφαρμογή σε σύνδεση ακριβείας/Принасовка насадки

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GB User Instructions Attachment bond®

Attachment bond is a composite luting resin which consists of a base paste and a catalyst. Curing is possible by light polymerization, with a color change as polymerization indicator, or by the chemically initiated autopolymerization. When autocured, the color change of Attachment bond will take some days.

Contents

Catalyst: Dimethacrylate, inorganic Fillers, Dibenzoylperoxide
Basepaste: Dimethacrylate, inorganic Fillers

Area of application

– Bonding of custom made attachment components to case frameworks (e.g. telescope crowns, tapered crowns).
– Bonding of off-the-shelf anchors, implant components and precision attachments to cast frameworks.

Design considerations for luting joints

The luting joints have to be designed strongly and with sufficient retention. The width of the luting gap should range between 0.1 and 0.2 mm only, which minimizes the luting stress and optimizes the bond strength.

Surface preparation and luting

When blocking out and preparing for duplication, make sure the joining surfaces lock circularly. The luting gap must not exceed 0.2 mm. When male and female parts fit passively on the model, the luting surfaces are air abraded with clean, new Al₂O₃ (cast frameworks with 250 µm and 4 bars pressure), attachments with 110 µm and 2 bars pressure). After air abrading do not touch the surfaces with anything. Using a Kulzer Resin-Metal-Bondingsystem is recommended (see the separate user instructions).
If a Resin-Metal-Bondingsystem is not available, clean the sandblasted surfaces with fresh acetone or ethanol and a clean brush. Protect with block-out wax all areas of the attachment which must not come in contact with the luting cement. Depending on the size of the luting surfaces, mix equal parts of base paste and catalyst avoiding air entrapment. Apply the cement mix to the luting surfaces with a spatula and avoid touching the surfaces with the spatula. Fix the attachment in its final position; total working time at normal room temperature is approximately 5 min

Polymerization

If light polymerization is used, align the luting joint to the light source for best exposure; light cure from all sides from which luting gap is accessible. The following exposure times are for each accessible side of the cement gap.

Curing unit	Polymerization time	Pre-cure polymerization time
Unilux® AC	BLZ 25	
Dentacolor® XS/UniXS®	90 s	
Translux EC/CL		20 s
Autopolymerization (23 °C)	20 min	

Finishing

Remove excess and finish after complete polymerization with standard carbide burrs.

Storage advice

Storage in the refrigerator (8 °C) is absolutely necessary. Before use, let Attachment bond warm-up to room temperature for approximately 30 min. Do not use material beyond its expiration date and keep away from sunshine. Turn the threaded plunger back a quarter revolution and put the caps back onto the syringes immediately after use. Store away from children's reach.

Safety advice

Irritating to eyes, respiratory system and skin. May cause sensitization by skin contact. Avoid contact with eyes. In case of contact with eyes rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water and soap.

Advice

The luting gap width must not exceed 0.2 mm. If the parts cannot be luted right after conditioning (30 min) the surfaces need to be sealed with a very thin layer of c & b liquid and light cured for 90 s in the Dentacolor XS unit, UniXS unit or Unilux AC BLZ 15. Attachment bond is stable up to 110 °C. Debonding of Attachment bond composite for repairs can be achieved by heat impact; electric spot welders are very useful for this procedure. Luting joints without a sufficient overlap and without mechanical interlock are not recommended.

Delivery units

Combination package
2 x 3 g net in syringes
Base- and catalyst paste

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® = Registered trademark of Heraeus Kulzer GmbH
Heraeus Kulzer GmbH