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1. Fields of application

Heat-curing resin for producing complete and partial dentures in the flask and injection technique. For permanent and removable dentures.

2. Safety precautions

The product contains polymerisable monomers (e.g. meth/acrylates) which can cause skin sensitisation in sensitive persons. In case of allergic reactions, stop contact with the product immediately. Do not use in case of known allergies to meth/acrylate-based resins.

Liquid:

Danger

H225 Highly flammable liquid and vapour.
H335 May cause respiratory irritation.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.

2.1 Prevention

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

2.2 Response

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P370+P378 In case of fire: use carbon dioxide (CO₂), sand or extinguishing powder to extinguish fire.

3. Processing instructions

3.1 Preparation

A hard modelling plaster type III is recommended for embedding the denture set up in wax, e.g. FINO MODEL STONE or FINOPLASTER. Carefully remove wax residues from the plastic components. Then insulate the wetted plaster model twice with a thin coating of suitable separation agent, i.e. FINOSOL alginate insulation (Item no. 11445). Observe the information provided in the instructions for use of the separating agent.

3.2 Processing

FINODON can be dosed freely, recommended ratio powder : liquid 10 g : 4 ml or multiples.

Add liquid to a mixing bowl. Angle and rotate the mixing bowl while sprinkling powder into the liquid until a non-wetted layer of powder remains. Tap the mixing bowl slightly when sprinkling the powder. Shake off excess powder after 30 s. Mix the dough well with a spatula and allow to swell for 15 mins. in the covered mixing bowl.

Optimal consistency is achieved when the material can be removed easily from the mixing bowl, when it no longer sticks and the mass no longer draws threads. Do not add any further powder. The processing time of FINODON depends on the room temperature and the mixing ratio and may vary for free dosing.

3.3 Packing/pressing process

Place resin dough together with excess in the flask, close flask with a press slowly at intervals and leave under pressure for approximately 10 mins. Add small amounts of dough for larger dentures. Recommended flask temperature is between room temperature and 45 °C.

3.4 Polymerisation

3.4.1 Short-term polymerisation

Ensure that a temperature of 80 °C is maintained for at least 15 mins. in the flask. Place flask in boiling water and boil for 20 mins. Then allow flask to cool slowly on the bracket. In case of thick denture parts, turn off the heat source after 20 mins. of boiling time and boil again for 20 mins. after a waiting period of 20 mins. Allow the flask to cool slowly in the water bath.

3.4.2 Water bath polymerisation

Place flask in at least 3 litres of lukewarm water. The flask must be completely covered with water. Heat water bath to 70 °C over approx. 30 minutes and maintain temperature for 30 mins. Then bring water bath to the boil over approx. 20 mins. and maintain boiling temperature for at least 30 mins. Allow the flask to cool slowly in the water bath.

3.4.3 Long-term polymerisation

Place flask in a heating bath filled with cold water and polymerise for 10 hours at 90 °C. Allow the flask to cool slowly in the water bath.

Note:
Keep dentures fabricated with FINODON in water until insertion.

3.5 Repairs

Execute repairs with FINODUR or FINOPRESS. Observe the instructions for use of the cold polymer.

4. Storage

P403+P233 Store in a well-ventilated place. Keep container tightly closed.

P410 Protect from sunlight.

P411 Store at temperatures not exceeding 21 °C.

P102 Keep out of reach of children.

Do not use beyond expiry date.

5. Physical data

Mixing ratio powder : liquid 10 g : 4 ml; or free dosing, colour powder and liquid: pink. Dough swell time 15 minutes, processing time depends on the room temperature and the mixing ratio and may vary for free dosing.

Polymerisation time:

Short-term polymerisation:

Boiling time 20 mins.,

In case of thicker sections:

Boiling time 2 x 20 mins.,

Steep time (waiting period between boiling times) 20 mins.,

Water bath polymerisation:

Heating time to 70 °C 30 mins.,

Hold time at 70 °C 30 mins.,

Heating time to boiling temperature 20 mins.,

Hold time boiling temperature at least 30 mins.,

temperature

Long-term polymerisation:

at 90 °C 10 h;

Cooling time slowly in water bath.

5.1 Chemical composition

Powder based on polymethyl methacrylate, liquid based on methyl methacrylate/dimethyl acrylate. Free of cadmium.

6. Delivery forms

FINODON denture resin		
powder, pink	1000 g	02412
liquid, pink	500 ml	02421

7. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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