poli
HOT NEVE

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PoliDent

poli
HOT NEW

POLIHOT is an acrylate for the fabrication of full and partial dentures by using the heat curing.

Colour range: Pink, Pink Veined.

Indications:

- full dentures
- partial dentures

Directions for use

- Preparation
After the complete washing out of the wax from the flask, the surface of model is to be isolated by an isolating fluid and thoroughly dried. In order to ensure a tight attachment to the base plate of the prosthesis apply the Polihot liquid on the roughened teeth.
- Mixing and packing/pressing
The recommended mixing ratio is 22 g powder: 10 g liquid (30 ml powder : 11 ml liquid).
Place the liquid into a mixing jar and add the powder little by little. The paste has to be thoroughly mixed and then put aside in a covered mixing jar for **3 minutes** at room temperature of 22°C. After that the paste has to be knead well and placed into a flask which is to be put into the pressure vessel.
WARNING: The paste has to be pressed into the flask within **7 minutes** after the liquid and powder have come into contact, otherwise the paste results in being too hard for pressing. The dough time as well as the working time are subject to the room temperature as well as the temperature of the liquid and the powder before mixing.
- Polymerisation
Place the flask in the mounting device into a water bath filled with cold water, heat the water to 100°C and polymerize for 35 - 45 minutes.
- Cooling
Cool the flask down to room temperature.
- Trimming
Open the cooled down flask and remove the stone. The prosthesis is trimmed and polished according to procedures common in dental laboratories.

- Repairs
The use of the cold curing material Policold is recommended, the colours of which match the acrylate for the fabrication of prosthesis by using the heat curing (Polihot).

Storage

Store at room temperature up to 25°C. Avoid direct exposure to sun.

Warning

Close container immediately after use. Do not use material after expiry date. Keep out of reach of children. Avoid direct contact between skin and uncured material as well as inhalation of the liquid.

Warning of liquid componente

Liquid is highly flammable. Irritating to eyes, respiratory system and skin. May cause sensitization by skin contact. Keep container in a well ventilated place. Keep away from sources of ignition. NO smoking. Do not empty into drains. Take precautionary measures against static discharges.

Packaging

02866 Polihot powder - 700 g, measuring cup
02867 Polihot liquid - 500 g, measuring cup

PoliDent