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English

**Instructions for use**

**Please read carefully before use**  
**Keep for referral purposes**

**1. Fields of application**

Universal cold-curing polymer for fabricating full and partial dentures as well as for repairs and linings. Fabrication using resin casting and injection techniques. For permanent and removable dentures for extraoral application. Direct lining is contraindicated.

**2. Safety precautions**

The product contains polymerisable monomers (e.g. meth/acrylates) which can cause skin sensitisation in sensitive persons. In case of allergic reactions, stop contact with the product immediately. Do not use in case of known allergies to meth/acrylate-based resins.

Liquid:

Danger

H225 Highly flammable liquid and vapour.

H335 May cause respiratory irritation.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

**2.1 Prevention**

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P241 Use explosion-proof electrical/ventilating/lighting/equipment.

**2.2 Response**

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P321 Specific treatment (see on this label).

**3. Processing instructions**

**3.1 Preparation**

A hard plaster type III is recommended for preparing the model, e.g. FINOPLASTER (Item no. 10101) or FINO MODEL STONE (Item no. 10251) hard plaster for modelling. Do not use any alcohol-based plaster-to-wax separating agents for the wax set-up as these could lead to white discolouration on the basal surface of the denture.

Secure the teeth set up in wax using a matrix made, for example, from FINOPASTE (Item no. 15180), S80 (Item no. 15191), S90 (Item no. 15192), DUO (Item no. 15175), DUO 95 (Item no. 15177) silicone putty or using a casting flask, for example, FINOFLOW casting flask (Item no. 14005). Then carefully remove wax residues from the plastic components. Subsequently insulate the plaster model carefully with a suitable separation agent, i.e. FINOSOL Alginate Insulation (Item no. 11445). Observe the information provided in the instructions for use of the separating agent. Remove excess separating agent carefully with compressed air. To improve the bond between FINOPRESS and acrylic teeth, roughen with a coarse diamond grinding instrument (remove grinding dust) and wet basal area with bonding agent, e.g. FINOCRYL BOND bonding agent (Item no. 02438).

**3.2 Processing**

**3.2.1 Casting technique**

FINOPRESS can be dosed freely, recommended ratio powder : liquid 10 g : 7 ml.

Add liquid to a mixing bowl and sprinkle quickly with powder. Mix immediately with a spatula for 15 s.

Further addition of powder in small amounts only and quickly. Shake mixing bowl to avoid formation of bubbles. FINOPRESS can be cast for approx. 3 min after mixing. After approx. 4 min start of the plastic phase, duration approx. 7 min.

**3.2.2 Polymerisation**

In the pressure pot, no earlier than 7 mins., no later than 13 mins. The surface of the resin should look dull.

Polymerisation temperature 55 °C, Polymerisation pressure 2 bar, Polymerisation time full denture 30 mins., all other work 20 mins.

**3.2.3 Injection technique**

FINOPRESS can be dosed freely, recommended ratio powder : liquid 30 g : 15 ml (2 : 1 corresponds to the amount for a full denture of average size). Add powder and liquid to a mixing bowl. The time point for injection is reached as soon as the mixed dough obtains a dull surface. The processing time depends on the ambient temperature. For fabrication of dentures using the injection technique, please observe the instructions for use provided by the manufacturer of the injection system.

**3.2.4 Polymerisation**

In a pressure pot:  
polymerisation temperature 55 °C,  
polymerisation pressure 2 bar,  
polymerisation time 30 min.

**3.3 Work up and polishing**

Following re-occlusion, remove the denture from the model and work-up with cross-linked carbide milling cutters. Smoothen the denture surface, e.g. with FINOCORUND grinding paper or flint paper. Polish the denture with pumice powder, e.g. FINOBIM. Then bring to a high polish with a high-gloss polishing agent suitable for prosthetic resins.

After manufacturing, all MMA/PMMA-based prosthetic materials contain a permissible residue of monomer which can lead to possible intolerance. To reduce the risk, place the denture made with FINOPRESS for a min. of 12 h in lukewarm water prior to insertion in the mouth.

**3.4 Extensions and repairs**

Execute extensions and repairs with FINOPRESS, FINODON or FINODUR.

**4. Storage**

P403+P233 Store in a well-ventilated place. Keep container tightly closed.

P410 Protect from sunlight.


P411 Store at temperatures not exceeding 21 °C.

P102 Keep out of reach of children.

Do not use beyond expiry date.

**5. Physical data**

Mixing ratio powder : liquid 10 g : 7 ml; or free dosing, colour powder and liquid: pink, mixing time casting technique max. 13 mins., time for casting up to 3 mins., in plastic state up to max. 7 mins., processing time for injection technique variable.

Important!

All data apply to a room temperature of 23 °C. Higher temperatures shorten the processing time, lower temperatures prolong the processing time.

**5.1 Chemical composition**

Powder based on methyl methacrylate copolymer

liquid based on methyl methacrylate/dimethylacrylate.

Cadmium-free.

The catalyst does not contain tertiary amines.

**6. Delivery forms**

FINOPRESS denture resin		
powder, pink	100 g	02434
	1000 g	02436
liquid, pink	80 ml	02435
	500 ml	02437

**7. Guarantee**

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

**8. Disposal**

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

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