
3M ESPE

Visio™ Beta vario

Light Unit for 3M ESPE Sinfony™ Polymerization

- ⓓ **Lichtgerät zur Polymerisation von 3M ESPE Sinfony™**
- ⓕ **Appareil servant à photopolymériser 3M ESPE Sinfony™**
- Ⓡ **Fotopolimerizzatore per 3M ESPE Sinfony™**
- ⓔ **Aparato de fotopolimerización para la polimerización de 3M ESPE Sinfony™**
- Ⓟ **Fotoemissor para a polimerização de 3M ESPE Sinfony™**
- ⓃⓁ **Polymerisatieapparaat voor 3M ESPE Sinfony™**
- ⓖⓇ **Συσκευή για τον φωτοπολυμερισμό του 3M ESPE Sinfony™**
- Ⓢ **Ljushärdningsapparat för polymerisation av 3M ESPE Sinfony™**
- ⓕⓂ **Valolaite 3M ESPE Sinfony™:n polymerointiin**
- ⓓⓀ **Polymeriseringslampe til 3M ESPE Sinfony™**
- Ⓝ **Lysapparat til polymerisering av 3M ESPE Sinfony™**

Operating Instructions

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Bedieningshandleiding

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1. Safety

PLEASE NOTE!

Prior to installation and start-up of this unit, carefully read the instructions provided herein! As with all technical devices, the proper function and safe operation of this unit depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations provided in these Operating Instructions.


1. The unit must be operated in strict accordance with the following instructions. The manufacturer accepts no liability for any damage resulting from the use of this unit for any purpose other than specified.
2. Prior to start-up of the unit make sure that the operating voltage stated on the rating plate is compatible with the available power supply. Operation of the unit at a different voltage may damage the unit.
3. Condensation resulting from the unit being transferred from a cold to a warm environment may be a potential hazard. Hence, the unit should be turned on only after it is completely equilibrated to ambient temperature.
4. Operate the device only with an appropriately rated power outlet with protective conductor. If you use an extension cord, please make sure that the protective conductor is not interrupted.
5. Prior to any change of lamps, the unit must be turned off and disconnected from the power supply.
6. In order to avoid electrical shock do not introduce any objects into the unit with the exception of replacement parts handled in accordance with the Operating Instructions.
7. Replace defective parts exclusively with original 3M ESPE replacement parts in strict compliance with the Operating Instructions. The product guarantee does not cover any damage resulting from the use of third-party replacement parts.
8. Should you have any reason to suspect the safety of the unit to be compromised, the unit must be taken out of operation and labeled accordingly to prevent third parties from inadvertently using this possibly defective unit.
The safety may be compromised, e.g. if the unit malfunctions or is visibly damaged.
9. Solvents, flammable liquids, and powerful heat sources must not be kept in the vicinity of the unit as either of these factors may damage the plastic housing of the unit.
10. Do not allow cleaning agents to enter the unit as this may short-circuit the unit or cause potentially dangerous malfunction.
11. Only factory-authorized technicians are authorized to open the housing of the unit and carry out any repairs.
12. It is not permissible
 - to cover the unit with a glass plate instead of the designated lid (may lead to bursting since there is no illumination from above),

Table of Contents	Page
1. Safety	1
2. Product Description	2
3. Technical Data	2
4. Start-up	2
5. Operation	2
6. Malfunction	3
7. Maintenance	3
8. Care	4
9. Storage / transport	4
10. Customer Information	4

- to plug or clog the air outlet opening inside the polymerization chamber (may produce a smear layer),
- to bypass the lid contact switch (may lead to insufficient illumination and vacuum),
- to place objects inside the polymerization chamber which are not specifically intended for this purpose (may damage the unit),
- to open the lid by force, e.g. with a screw-driver (may damage the unit).

2. Product Description

Visio Beta vario, manufactured by 3M ESPE, is an illumination unit for the polymerization of composite veneer

materials and light-curable dental materials, such as light-curable tray materials, modeling resins for custom-made attachment parts or block-out materials. There are four programs with different vacuum periods available which can be selected to suit the specific polymerization requirements of the material in use. In the final step of the procedure, in which the material is exposed to a vacuum, the oxygen inhibition layer commonly observed with composite materials is eliminated.

For details on all products mentioned below please refer to the corresponding Instructions for Use.

Operating Instructions should not be discarded for the duration of product use.

3. Technical data

Supply voltage:	120 V, 230 V \pm 10% 50/60 Hz pre-set voltage and AC-frequency are stated on the rating plate
Power consumption:	approx. 250 W (incl. vacuum pump)
Dimensions:	280 x 360 x 135 mm (D x W x H)
Weight:	approx. 6.4 kg
Ambient temperature range:	18°C/64°F... 40°C/104°F
Relative humidity of air:	max. 80% at 31°C/88°F max. 50% at 40°C/104°F

Operating the unit with a cold pump may lead to overloading of the pump and damage the pump due to poor lubrication!

For details on the installation of the Visio Beta Vacuum pump, manufactured for 3M ESPE, please refer to separate operating instructions!

Push the suction hose of the Visio Beta Vacuum pump (labeled "Gerät-unit") onto the hose coupling piece of the Visio Beta vario unit.

Plug the power cord of the vacuum pump into the upper socket on the back of the Visio Beta vario unit. Connect one plug of the power cord to the lower socket on the back of the Visio Beta vario unit, and the other plug of the cord to a properly grounded power outlet.

4. Start-up of the unit

ATTENTION

Do not start-up the unit plus vacuum pump before the pump has equilibrated at least to room temperature.

5. Operation

Turn on the power switch. For a total of 3 sec the display shows the number of times the lamps have been activated (two digits show the lamps acts in units of thousands/hundreds) and both lamp error indication lights are on.

If the two lamp error indication lights flash and an acoustical signal (three sounds) is emitted, the lamps have been activated for a total of 9999 times (display shows 99) and the lamps must be replaced (please refer to the section on Malfunctions – “End of lamp life” and “Maintenance and Care”).

Subsequently, the display shows “Pu” (for a period of three seconds) and the red vacuum error indication lamp lights up. If you press the “Start” button during this time interval, the vacuum pump regeneration program (duration of 99 minutes) is initiated (for details please refer to the Visio Beta vacuum pump’s operating instructions). The unit is ready for use once all lamps are off.

Filling the polymerization chamber

The lid of the polymerization chamber opens once the “Lid” button is pressed. Materials to be polymerized may be placed in the polymerization chamber in any arrangement or order. Please ensure that the parts do not overlap! Then close the lid by folding it down.

ATTENTION

Crowns and bridges must be placed in the chamber without the plaster cast as the residual water content of the plaster prevents the required vacuum from being established, and may damage the vacuum pump. (Exception: part 1 of program 4).

Program selection

Use the buttons, “P1” through “P4”, to select one of the following four programs. The selected program is indicated by a green illuminated area:

- Program 1: 1 min illumination, 14 min illumination and vacuum
Application: Final polymerization of Sinfony, manufactured by 3M ESPE
- Program 2: 7 min illumination, 10 sec illumination and vacuum
Application: Opaquer polymerization for multiple veneer units or for coating cast retentions with Opaquer Pink
- Program 3: 1 min illumination, 1 min illumination and vacuum
Application: Polymerization of cold light-curable dental materials at a wavelength between 400 and 500 nm
- Program 4: 4 min illumination, pause, 4 min illumination and vacuum
Application: Polymerization of light-curing tray materials During the break (display shows “2.-”, three short beeps are emitted) open the lid and remove the plaster cast. The second part of the program is initiated by pressing the “Start” button.

Start-up of the program

After selecting the desired program, press the “Start” button upon which the program proceeds automatically. (If you fail to start up the program within 10 seconds

after selecting any of programs 2, 3 or 4, the device automatically returns to program 1 for safety reasons.)

During the execution of the program the remaining time is displayed in units of minutes (in units of seconds in the second part of program 2). All programs can be terminated by pressing the “Stopp” button.

After completion of the program, the lid of the polymerization chamber opens automatically, the signal indicating the end of the program (5 short beeps) is emitted, and the unit automatically returns to program 1. The work pieces may now be removed from the chamber and processed without delay.

6. Malfunction

Both the vacuum and the lamps are controlled by the unit. The user is notified of any operating errors by means of the three red error indication lights and an acoustical signal.

The vacuum pump shuts itself off if the required vacuum fails to be established during polymerization. The unit should not be shut off at this point, since the polymerization process can be continued for the remaining polymerization time (by pressing the “Start” button) once the source of the error has been eliminated.

Controls to be carried out in case of a vacuum error:

- Check the hose between the unit and the vacuum pump for leaks
- Check the rubber gasket for leaks and the glass cover plate of the polymerization chamber for contamination.

Controls to be carried out in case of a lamp error:

- Restart the program to make sure the operating error is not caused by a single short failure of the lamp to light up. If the operating error persists, the pair of lamps must be replaced (please refer to “Maintenance and Care”).

End of lamp life:

The two red lamp error indication lights and a triple acoustical signal are activated once the lamp has been turned on for a total of 9999 times (display shows 99). At this point the lamps need to be replaced (please refer to “Maintenance and Care”) since the performance of the lamps may no longer be sufficient to ensure proper polymerization.

This signal is emitted each time the power button is activated and each time a program ends or is terminated.

After replacement of all four lamps, the lamp activation counter must be reset to zero. In order to do so, turn the unit off, then press the buttons, “Start/Stop” and “Lid” (keep depressed), and turn on the power switch.

7. Maintenance

ATTENTION


In order to avoid any wear-related loss of lamp performance and to ensure proper polymerization, all four Visio Beta fluorescent lamps LB 1, manufactured for 3M ESPE,

must be replaced after 9999 illumination cycles (see "Malfunction").

Replacement of the upper pair of lamps
(inside the cover of the polymerization chamber)

Open the polymerization chamber by pressing the "Lid" button.

Turn the unit off and pull the power plug from its socket!

Unscrew the six screws located in the cover plate, and remove both the cover plate and glass plate.

Clean both sides of the glass plate with a cloth dampened in alcohol.

Remove the two lamps from their respective holders by pulling in an

upward direction. If the reflector located behind the lamps is contaminated, clean it with a cloth dampened in alcohol.

Caution: Do not scratch the reflector with any hard objects!

Place new fluorescent lamps LB1 in the respective holders.


Replace the cover plate and fix in place by tightening the six screws.

Replacement of the lower pair of lamps


Empty the polymerization chamber and close the lid.

Turn the unit off and pull the power plug from its socket!

Place the unit on its right side.

Unscrew the screw located in the bottom door and remove the door.

Remove the two fluorescent lamps from their respective holders by pulling in an upward direction.

Clean the bottom door reflector and the bottom side of the glass cup with a cloth dampened in alcohol.

Caution: Do not scratch the reflector with any hard objects!

Place new fluorescent lamps LB1 in the respective holders.

Attach bottom door by its hinges, close the door and tighten the screw.

8. Care of the unit

Clean the unit with a commercial dish-washing detergent. (Organic) Solvents may damage the housing of the unit.

Any contamination inside the polymerization chamber should be removed with alcohol. Persistent contaminating material should be soaked by placing an alcohol-saturated piece of cloth inside the chamber and closing the cover. (Do not press the Start-Stopp button.) After extended soaking (preferably over night) use a plastic spatula to remove residues of the working materials from the polymerization chamber and glass plate. Do not use a metal instrument for this purpose!

Fluorescent lamps must be cleaned/dusted at least at semi-annual intervals to ensure proper performance of the lamps. For details on how to clean reflectors and fluorescent lamps, please refer to the section on "Replacement of upper/lower pair of lamps".

9. Storage / transport

Do not store Visio Beta vario below $-10^{\circ}\text{C}/24^{\circ}\text{F}$ or above $+60^{\circ}\text{C}/140^{\circ}\text{F}$.

10. Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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3M, ESPE, Sinfony and Visio are trademarks of 3M or 3M ESPE AG.
© 3M 2003. All rights reserved.
3M ESPE AG
3M ESPE Technical Hotline/MSDS Information in U.S.A. and Canada 1-800-634-2249.