FINOFRAME PMMA SPLINT

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- | | |
|-------------------------------|-----------------------------|
| (en) Instructions for use | (lt) Naudojimo instrukcija |
| (fr) Notice d'utilisation | (lv) Lietošanas informācija |
| (it) Istruzioni per l'uso | (et) Kasutusjuhend |
| (es) Instrucciones de uso | (hu) Használati utasítás |
| (pt) Instruções de utilização | (cs) Návod k použití |
| (pl) Instrukcja użycia | (sl) Navodila za uporabo |
| (nl) Gebruiksaanwijzing | (sk) Návod na používanie |
| (sv) Bruksanvisning | (bg) Инструкции за употреба |
| (da) Brugsanvisningen | (el) Οδηγίες χρήσης |
| (no) Bruksanvisningen | (tr) Kullanım kılavuzu |
| (fi) Käyttöohje | |

FINOFRAME PMMA SPLINT

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Stand der Information 2019-10

FINOFRAME PMMA SPLINT

Milling blank made of thermoelastic PMMA without fillers and peroxide

en Instructions for use, please read carefully!

Intended use

- FINOFRAME PMMA SPLINT is suitable for the fabrication of occlusal splints, drilling templates, functional impression trays and impression trays, bite templates, bleaching trays and "lost moulds" for casting in dental prosthetics.

Indication

- Abnormal jaw function
- Total or partial tooth loss

Contraindication

- Hypersensitivity to one of the components

Composition

- Polymethylmethacrylate (PMMA)
- Methylmethacrylate (MMA) may be contained as residual monomer up to max. 1 %.

Physical data

Property	Result	Test method
Glass transition	121 °C	ISO 11357
Flexural strength	156 ± 6 MPa	EN ISO 10477
Modulus of elasticit	3000 ± 100 MPa	EN ISO 20795-1
Vickers hardness	23.2 ± 0.3 HV 0.2	EN ISO 6501
Water uptake	< 32 µg / mm ³	EN ISO 20795-1
Water solubility	< 1.6 µg / mm ³	EN ISO 20795-1

CAM processing

- In dry processing in the milling machine, care must be taken that milling chips are extracted or evacuated.
- FINOFRAME PMMA SPLINT can generally be used with the common burs available for thermoplastics and PMMA. If you have questions on these topics, please contact the manufacturer of your milling machine, software or milling tools.
- FINOFRAME PMMA SPLINT features a glass-transition temperature of 121 °C. Particularly in the case of dry milling, it should be ensured that there is no excessive heat transfer to the material, e.g. due to excess speed and/or feed rate, as otherwise the milling tool and/or the material surface could smear.

Trimming

- Cross-cut carbide burs suitable for resin should be used to rework and separate frameworks, crowns, bridges, denture bases, bite splints and drilling templates from the blank.
- To prevent inaccurate fitting after fabricating the construction, avoid excessive heat development during trimming and polishing.

Restoration

Constructions, which have been fabricated from the FINOFRAME PMMA SPLINT can be restored and supplemented at any time with all standard MMA-based cold-curing polymers. Prior conditioning of the surface is necessary, preferably with acrylate-based bonding agent in accordance with the specifications of the manufacturer.

Storage

- No special measures required.

Notes

- Minimum wall thickness after polishing ≥ 1 mm.
- Record the LOT number with every process that requires identification of the material.

- To ensure traceability and assignment of the product data as well as the LOT numbers for the patient at any time, it is recommended to store the milling blank in the original packaging during two milling operations.
- When disinfecting the construction, only alcohol-free disinfectants should be used as crazing could otherwise form.
- Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.

Warnings

Processing PMMA frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from **www.fino.com**

Disposal

Residues from processing the FINOFRAME PMMA SPLINT have to be disposed of as residual waste.

The product properties are based on compliance with and observation of these instructions for use.

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FINO GmbH
Merz Dental GmbH