

11

ALL-CERAMIC DENTIN SHADED INGOTS • LINGOT TEINTÉ POUR LE TOUT CÉRAMIQUE • LINGOTTI TONALITA' DENTINA PER CERAMICA INTEGRALE • ALL-CERAMIC DENTINFARBIGE PELLETS • LINGOTE COLOREADO DE DENTINA COMPLETAMENTE DE CERÁMICA

12

ALL-CERAMIC VALUE SERIES SHADED INGOTS • LINGOTS TEINTÉS VALUE SERIES POUR LE TOUT CÉRAMIQUE • LINGOTTI TONALITA' SERIE VALUE PER CERAMICA INTEGRALE • ALL-CERAMIC VALUE SERIES FARBIGE PELLETS • LINGOTES COLOREADOS VALUE SERIES COMPLETAMENTE DE CERÁMICA

MANUFACTURER
DENTSPLY International
DENTSPLY Ceramco
570 West College Avenue
York, PA 17401 U.S.A.
1-800-243-1942
prosthetics.dentsply.com

EC REP DeguDent GmbH
Rodenbacher Chaussee 4
63457 Hanau-Wolfgang
Germany
+49/6181/59-50


Distributed by/Distribué par/Distribuido da/
Vertrieb durch/Distribuido por:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3
905-851-6060

Directions for Use
Mode d'emploi
Istruzioni per l'uso
Gebrauchsanleitung
Instrucciones de Uso


1

PASTE OPAQUE APPLICATION • APPLICATION DE L'OPAQUE EN PÂTE • APPLICAZIONE DI OPAQUE IN PASTA • APPLIKATION DES PASTENOPAKERS • APLICACIÓN DE OPAQUE EN PASTA


2

PFM MARGIN PORCELAIN APPLICATIONS • APPLICATION DE LA PORCELAIN D'ÉPAULEMENT PFM • APPLICAZIONE DI DELLA MASSA MARGINAL CON LA TECNICA DELLA CERAMICA COTTA SU METALLO • APPLIKATION DER PFM-KERAMIKSCHULTER-MASSE • APLICACIÓN DE PORCELAINA PARA MÁRGENES PFM

3

BASIC BUILD-UP • MONTAGE ÉLÉMENTAIRE • STRATIFICAZIONE BASE • EINFACHER AUFBAU • CONSTRUCCIÓN BÁSICA


4

ADVANCED BUILD-UP • MONTAGE VANCÉ • STRATIFICAZIONE AVANZATA • FORTGESCHRITTENER AUFBAU • CONSTRUCCIÓN AVANZADA


5

CHARACTERIZATION OPTIONS • OPTION DE CARACTÉRISATION • OPZIONI DI CARATTERIZZAZIONE • CHARAKTERISIERUNG-SOPTION • OPCIÓN DE CARACTERIZACIÓN


6

PRESSED & STAINED VENEERS • FACETTES PROTHÉTIQUES PRESSÉES ET COLOREES • FACETTE PRESSATE E PIGMENTATE • GEPRESSTE UND EINGEFÄRBT VERBLENDUNGEN • CARILLAS TINTADAS Y PRENSADAS


7

FELDSPATHIC VENEERS • FACETTES PROTHÉTIQUES FELDSPATHIQUES • FACETTE REALIZZATE SU REFRATTARIO • FELDSPATVERBLENDUNGEN • CARILLAS FELDESPÁTICAS

8

PTM MARGIN • PTM ÉPAULEMENT • MARGINE CON LA TECNICA DELLA CERAMICA PRESSATA SU METALLO • PTM SCHULTERMASSE • PTM MARGEN

9

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Porcelain Application Directions for Use

INDICATIONS FOR USE

The Ceramco iC Porcelain system is indicated for:

1. Porcelain Fused-to-Metal Applications: Single and multiple unit porcelain fused-to-metal fixed prosthodontic restorations.
2. Press-to-Metal Applications: Single and multiple unit press-to-metal fixed prosthodontic restorations.
3. All-Ceramic Applications: Single unit anterior and posterior premolar metal-free fixed prosthodontic restorations, laminate veneers, inlays and onlays.

CONTRAINDICATIONS

Only the indications listed above are suitable.

WARNINGS

Consult applicable DENTSPLY SDS (Safety Data Sheets) for safe handling and usage. Use only with recommended Ceramco iC system components. Use protective equipment to guard against dust when grinding.

PRECAUTIONS

Refer to the special notes throughout the instructions for best practices for using Ceramco iC materials. Other porcelains are not thermally compatible and will not have the correct dilatometric characteristics, and therefore should not be used. While initial results with some materials may appear acceptable, internal stress can compromise long term success.

ADVERSE REACTIONS

The Ceramco iC system contains components which may cause skin dryness, sensitization (allergic contact dermatitis or other allergic reactions in susceptible persons).

Symbols on product labels:

Manufacturer

Batch Code

Reorder Number

Manufacture Date

Expiration date

Consult Directions For Use

Caution

Keep dry

STORAGE

Store materials in packages with the lid tightly sealed. Do not leave open bottles where dust and debris may contaminate the porcelain. Replace cap when not in use. Avoid vibration and heat. Keep dry. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence. Do not use after expiration date.

ALLOY SELECTION

Ceramco iC Porcelain is compatible with high noble, noble and predominately base alloys. Be sure to consult the alloy manufacturer for alloy compositions and coefficient of thermal expansion data. Ceramco iC will resist discoloration on silver containing alloys; however, routine purging of your porcelain furnace and sagger trays is highly recommended. Alloys should have a CTE of 13.8 – 15.1 @ 500°C, Ideal range: 14.0 – 14.6 @ 500°C.

LIQUIDS

Caution should be used in selecting modeling and staining liquids. Ceramco iC liquids are recommended for their ability to improve the handling and carving of porcelain, to reduce slumping during build-up, and prevent any silver discoloration. Distilled water may also be used. However, other build-up and/or stain liquids may leave a residue behind that could cause discoloration.

INSTRUCTIONS

PASTE OPAQUE APPLICATION (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

Note: For Compatis[®] CoCr the use of NE-Bonder is required prior to applying the paste opaque. For other nonprecious frameworks NE Bonder is recommended but not obligatory.

1. Oxidize or degass and clean the alloy according to the alloy manufacturer's instructions.
2. Rinse a flat brush in distilled water and wipe it on a dry paper towel. Check all castings to make sure that all surfaces are dry.
3. Dispense a small amount of the paste opaque onto a pallet. Water must not be mixed into the paste. If necessary, the paste may be thinned with Paste Modifier Liquid. Use a flat brush and using a light touch apply an even coat of the paste onto the metal substructure. Tap the hemostat lightly to even out the brush strokes. Achieve an even coat and do not allow the material to puddle. Opaque modifiers would be applied at this point and blended into the paste opaque prior to the crystal application.

Section #1, Image #1: Apply the opaque completely covering the coping

Section #1, Image #2: Use a light vibration to smooth the brush marks

Section #1, Image #9: Blend the opaque modifier with the dentin shaded opaque

Section #1, Image #10: Apply the mixture to the coping

4. Use a hemostat and hold the restoration over the collection bowl. Sprinkle the crystals onto the paste. Tap once and do not vibrate. Use a dry brush to remove any crystals from the inside of the coping and any non-porcelain bearing areas. Use a brush at a 45° angle to remove any crystals that may remain in the margin area.

Section #1, Image #3: Sprinkle crystals on the surface

Section #1, Image #11: Apply the crystals and fire

5. Dry and fire according to the recommended temperatures. The correct visual indicator is a slightly shiny surface with a small amount of texture.

Section #1, Image #4: Fire according to the recommended schedule

Section #1, Image #12: Opaque and Opaque Modifier

Note: Drying too rapidly results in small voids in the paste under the crystals. These small voids may come to the surface during subsequent firings. Following the recommended dry and pre-dry times will ensure this does not occur. It is also effective to place the restoration on a sagger tray and put the assembly on a warm surface or hot plate prior to placing in the porcelain furnace.

Note: In most cases a single layer of paste opaque and crystals will provide sufficient coverage to mask out the alloy. However, if necessary a two coat application may also be utilized. In this case, apply the first coat of opaque BUT DO NOT APPLY ANY CRYSTALS. Fire the opaque, and follow with a second coat of opaque and crystals.

Section #1, Image #5: Apply the first coat of opaque and fire with NO crystals

Section #1, Image #6: Apply a second coat of opaque

Section #1, Image #7: Apply the crystals

Section #1, Image #8: Fire according to the recommended schedule

MARGIN PORCELAIN APPLICATION (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

SUBSTRUCTURE DESIGN

Extend the facial surface of the metal substructure so

that the metal ends at the axial gingival corner of the margin preparation or extends slightly onto the gingival ledge.

Section #2, Image #1: Alloy cut-back

Section #2, Image #2: Opaque applied and fire

DIE SEALER APPLICATION

1. Using the DENTSPLY Prosthetics Die Sealer, apply a thin, even layer of sealer over the die margin area. Blow off any excess die sealer and allow the liquid to thoroughly dry.
2. Reapply the die sealer at least 2 more times and allow to dry.
3. Mark the facial margin with a wax pencil.

DIE RELEASE APPLICATION

1. Using the DENTSPLY Prosthetics Die Release, apply a thin, even layer of die release over the margin area. Blow off the excess die release.
2. Apply and allow to dry several more times until the die has a very slight sheen.

FIRST MARGIN APPLICATION

1. Place a small amount of the shaded margin porcelain on a clean palette. If modifiers are necessary, add the modifier porcelain to the shaded margin porcelain.
2. Mix the Ceramco iC Margin and Opaque Corrector Liquid with the porcelain to a thick creamy consistency. This mixture will air harden in approximately 5 minutes. If additional working time is required, mix in a small amount of distilled water.
3. Place the porcelain mixture in the gingival area, pushing the mixture down to but not over the gingival margin. Do not overbuild.

Section #2, Image #3: First margin application

Note: All porcelains shrink by volume. The larger the porcelain mass the more apparent the firing shrinkage. Better results are achieved by applying a small amount of porcelain on the first porcelain application.

4. Allow the margin porcelain to dry completely (at least 5 minutes). The porcelain surface will be chalky and hard.
5. Carefully lift the coping from the die.
6. Dry and fire the restoration according to the recommended temperatures. The correct visual indicator will be shiny yet grainy.

Section #2, Image #4: First margin fired

SECOND MARGIN APPLICATION

1. Be sure the die is clean and reapply the die release. Allow to dry.
2. Place the coping on the die. Make a new mix of margin porcelain and margin liquid. Apply the mixture to the gingival margin.
Section #2, Image #5: Second margin application
3. Allow the margin porcelain to dry completely (at least 5 minutes) and remove the coping from the die.
4. Dry and fire the restoration according to the recommended temperatures.
5. The correct visual indicator will be shiny yet grainy.

Section #2, Image #6: Second margin fired

Note: If additional applications are required, repeat the procedures used for the second margin application.

CERAMCO IC BASIC AND ADVANCED BUILD-UP (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

OPACEOUS DENTINS

Natural teeth exhibit numerous color zones. These areas have various degrees of translucency, chroma and opacity. The experienced technician learns to use a variety of shades and modifiers to effectively mimic these characteristics. Frequently, the task is compromised by circumstances that do not allow for adequate porcelain thickness. The Ceramco iC Opaque Dentin porcelains address the problem of limited space, simplify internal color modifications and aid overall. The Ceramco iC Opaque Dentin porcelains are the same hue as the Ceramco iC Dentin shades, but have about 10% more opacity. They may be used undiluted, or mixed with the Ceramco iC Dentins or Dentin modifiers.

THIN AREAS

Frequently insufficient dentin/incisal thickness exists (less than 1.0 mm) over the opaque porcelain. The shade in this thin area can be improved by applying a layer of the opaque dentin porcelain under the same dentin porcelain shade. When contouring, you will be able to reduce the area where the opaque dentin was layered and the prescribed shade will be maintained. Complete the building and firing of the restoration in the normal manner.

Section #5, Image #1: Thin facial areas

FACIAL MARGIN (GINGIVAL COLOR)

If the metal substructure has been designed with a feather margin, the opaque porcelain around the facial margin leaves minimal room for the dentin porcelain, and the opaque porcelain will be highlighted. Esthetics will be improved by applying opaque dentin porcelain alone or as a blend with 10% to 20% of one of the dentin

modifiers. This technique creates a chromatic cervical third, common in natural dentition.

Section #5, Image #2: Gingival

PONTICS

A color imbalance frequently occurs between the pontics and the abutments of a multiple-unit fixed partial denture. The color difference is caused by the absence of a full metal substructure in the pontic area. Ceramco iC Opaque Dentin porcelain applied to the bottom and the cervical area of the pontic creates a uniform color completely across the fixed partial denture.

Section #5, Image #4: Bottom of pontics

LINGUAL FOSSAS AND OCCLUSAL SURFACES

Often, inadequate space is available on the lingual surface of an anterior unit, or on the occlusal surface of a posterior unit. Ceramco iC Opaque Dentin porcelain applied alone or as a blend with 10% to 20% dentin

modifiers in these thin areas prevents the opaque porcelain from being highlighted.

Section #5, Image #3: Lingual

CERAMCO IC HIGH CHROMA (HC) AND DENTIN PORCELAINS

The Ceramco iC High Chroma porcelains (HC) are dentin modifiers that have the same hue, value and opacity as the Ceramco iC Dentin porcelain but approximately 70% higher in chroma (color saturation). The HC porcelains can be used undiluted or mixed with the Ceramco iC Dentins, Opaque Dentins, Natural Enamels, or Opal Enamels for a variety of effects. For a very subtle effect, mix in a very small amount of dentin modifier with the dentin porcelain.

CERAMCO IC NATURAL ENAMELS, OPAL ENAMELS, AND MAMELON PORCELAIN

The Ceramco iC Natural Enamels are non-opalescent and can be used in thin or thick layers over the dentin to duplicate the vitality of natural teeth. They can also be mixed with the dentin or opal enamel porcelains for a variety of effects. The Ceramco iC Opal Enamels are an intense opal enamel porcelain. They can be mixed in with the natural enamels or used as an effect powder. The mamelon porcelains can be used in full concentration or diluted with the clear enamel porcelain. They are normally placed on the dentin porcelain lobes prior to the enamel porcelain application.

DENTIN BUILD-UP TECHNIQUE

1. Areas on the model that the porcelain will contact should be completely sealed with DENTSPLY Prosthetics Die Sealer. Other sealers or oils may leave a residue that will cause discoloration and should not be used.

- Mix the porcelains with Ceramco iC Modeling Liquid U or distilled water to a paste-like consistency.
- Build dentin porcelain or opaque dentin/dentin modifier mixture around the cervical one third. Condense slightly. If building a multiple unit fixed partial denture, apply the dentin porcelain mixture to the saddle area of the pontic and seat the framework on the model. Continue to apply the dentin porcelain until full contour has been achieved. Condense slightly.

Section #3, Image #1: Full contour Build-Up

Section #3, Image #2: Full contour Build-Up

Section #4, Image #1: High Chroma gingival, Opaque Dentin at the Incisal

- Cutback the incisal two-thirds to allow room for the enamel porcelains. If required, cutback to accentuate the lobe formation found in natural teeth. Be sure to remove adequate porcelain in the interproximal area.

Section #3, Image #3: Reduce the incisal 2/3rds down the facial. Maintain the incisal height.

Section #3, Image #4: Incisal Cut-Back

Section #3, Image #5: Incisal Cut-Back

Section #3, Image #6: Incisal Cut-Back

Section #4, Image #2: Dentin build-Up with Incisal Cut-Back

Section #4, Image #3: Incisal Cut-Back

- If mamelons are required, place a small amount of mamelon porcelain on the dentin lobes. To simulate natural dentition, vary the length of the mamelons.

Section #4, Image #4: Mamelon Porcelain

- Apply the enamel porcelain to the incisal or occlusal one-third, overbuilding by approximately 10%. If desired layer the selected enamel porcelains in narrow vertical columns across the facial-incisal area to simulate the enamel rod formations found in natural teeth. Use natural enamel porcelains of different colors, or the clear enamel porcelain, or the opal enamel porcelains arranged in discrete parallel columns to create contrast. Condense lightly.

Section #3, Image #7: Enamel porcelain application

Section #3, Image #8: Over build by 10%

Section #3, Image #9: Enamel application complete

Section #3, Image #10: View from the facial

Section #4, Image #5: Incisal and clear lobes

Section #4, Image #6: Super Clear on the facial

- Remove the coping from the model. Apply a small amount of enamel porcelain to the mesial/distal contacts. Slightly condense the porcelain.

Section #3, Image #11: Add enamel porcelain to the contacts and fire

- If building a multiple unit fixed partial denture, use a thin bladed instrument and cut slightly into the interproximals. It is not necessary to cut all the way down to the opaque.

- Dry the build-up away from direct heat for at least five minutes before bringing it close to the open furnace entrance. Increase the drying time to ten minutes for large multiple unit fixed partial dentures.

- Fire according to the recommended temperatures. A proper visual indicator is shiny with a small amount of surface texture.

Section #3, Image #12: First Bake – Shiny with some texture

Section #3, Image #13: Completed Crowns

Section #4, Image #7: Visual Indicator after firing is shiny with a little texture

Section #4, Image #8: Completed Crown

CERAMCO IC CHARACTERIZATION OPTIONS (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

Ceramco iC offers a number of materials and modifiers to create an infinite number of unique shading possibilities. All the components can be used interchangeably for PFM's, all-ceramic and press-to-metal restorations. Please refer to the inside cover of this instruction for some examples.

Section #5, Image #5: Mamelon and Clear enamel effect

Section #5, Image #6: Opal enamel and High-Chroma dentins

Section #5, Image #7: Bleached shade with standard enamel

Section #5, Image #8: Bleached shade with opal enamel and Super Clear enamel

Section #5, Image #9: Internal crack lines

Section #5, Image #10: Super Clear enamel, Opal enamel, Stained crack line and stained margin area

Section #5, Image #11: Mamelons, Super Clear enamel

Section #5, Image #12: Opaque modifier on the occlusal table

Section #5, Image #13: Pit and Fissure Staining

CONTOURING AND CLEANING

1. After firing, use fine-grit, non-contaminating stones, discs, or diamond burs to refine the anatomy. If no porcelain additions are needed, the case may be glazed after thorough cleaning.
2. If additions are required or before glazing, lightly blast the porcelain surface with unrecycled 50-micron aluminum oxide at 20 psi. Clean the case with a steam cleaner or ultrasonic bath using distilled water.

SECOND APPLICATION OF DENTIN AND ENAMEL PORCELAINS

1. Apply the additional porcelain using the same steps and techniques described for the initial application. Match the layering of the opal enamel porcelain over the dentin porcelain. Second applications of Ceramco iC porcelains require the same care in drying and firing as the first layer.

OPAQUE CORRECTION PORCELAIN

If during the contouring procedures the porcelain surface is removed and the metal surface is exposed, this may be easily repaired by using the opaque correction material. This correction opaque should be used for small minor repairs only.

1. Lightly blast the surface that is to be repaired with aluminum oxide at 20 psi. Then steam clean or clean in an ultrasonic cleaner in distilled water for 5 minutes.
2. Place a small amount of the correction opaque onto a pallet. Mix with a small amount of Ceramco iC Margin and Opaque Corrector Liquid to a paste-like consistency.
3. Apply the opaque mixture to the exposed metal and allow to dry for about 5 minutes. The opaque will appear chalky and hard.
4. Apply either Ceramco iC Opaqueous Dentin, Dentin or Add-On porcelain and fire according to the recommended firing temperatures.

ADD-ON PORCELAIN

The Ceramco iC Add-On porcelain is available in a translucent enamel, dentin light, medium and dark, and tissue colored pink, reddish-pink, salmon, and dark. They fire at slightly lower temperatures than the regular dentin porcelain. These porcelains are always fired under vacuum and are used when a low temperature add-on is required.

1. Lightly blast the surface that is to be repaired with aluminum oxide at 20 psi. Then steam clean or clean in an ultrasonic cleaner in distilled water for 5 minutes.
2. Mix the porcelains with Ceramco iC Modeling Liquid U or distilled water to a paste-like consistency.

3. Apply the mixture to the desired area.
4. Fire according to the recommended firing temperatures.

STAINS AND OVERGLAZE

1. If a natural glaze is required, be sure the porcelain surface is clean and follow the firing instructions for a natural glaze.
2. If an overglaze is required, mix the Ceramco iC Overglaze with the Ceramco iC Stain and Glaze Liquid and apply a thin layer of the mixture to the porcelain surface.
3. Apply the stain to the desired area.
4. Check the shade and add any necessary stains. Fire according to recommended temperatures.

FINAL MARGIN

Final Margin is used to make any minor corrections on a porcelain margin. They fire at a low temperature and should be used with the stain and glaze steps.

1. Place a small amount of the final margin porcelain on a clean palette.
2. Mix the Ceramco iC Margin and Opaque Corrector Liquid with the porcelain to a thick creamy consistency. This mixture will air harden in approximately 5 minutes. If additional working time is required, mix in a small amount of distilled water.
3. Place the porcelain mixture in the required area pushing the mixture down to but not over the gingival margin. Do not over build.
4. Allow the margin porcelain to dry completely (at least 5 minutes). The porcelain surface will be chalky and hard.
5. Carefully lift the coping from the die.
6. Dry and fire the restoration according to the recommended temperatures.

Laminate Veneers Directions for Use

INDICATIONS FOR USE

The Ceramco iC Porcelain system is indicated for:

1. Porcelain Fused-to-Metal Applications: Single and multiple unit porcelain fused-to-metal fixed prosthodontic restorations.
2. Press-to-Metal Applications: Single and multiple unit press-to-metal fixed prosthodontic restorations.
3. All-Ceramic Applications: Single unit anterior and posterior premolar metal-free fixed prosthodontic restorations, laminate veneers, inlays and onlays.

CONTRAINDICATIONS

Only the indications listed above are suitable.

WARNINGS

Consult applicable DENTSPLY SDS (Safety Data Sheets) for safe handling and usage. Use only with recommended Ceramco iC system components. Use protective equipment to guard against dust when grinding.

PRECAUTIONS

Refer to the special notes throughout the instructions for best practices for using Ceramco iC materials. Other porcelains are not thermally compatible and will not have the correct dilatometric characteristics, and therefore should not be used. While initial results with some materials may appear acceptable, internal stress can compromise long term success.

ADVERSE REACTIONS

The Ceramco iC system contains components which may cause skin dryness, sensitization (allergic contact dermatitis or other allergic reactions in susceptible persons).

STORAGE

Store materials in packages with the lid tightly sealed. Do not leave open bottles where dust and debris may contaminate the porcelain. Replace cap when not in use. Avoid vibration and heat. Keep dry. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence. Do not use after expiration date.

INGOT TYPES

The Ceramco iC system offers three different types of shaded ingots.

1. Dentin Shaded ingots come in all Vita Lumin[®], 3D-Master[®]*, and bleached shades. They are slightly higher in opacity and generally used with an enamel add-on technique for either press-to-metal, or all-ceramic applications.
2. Value Series ingots have a higher translucency, come in 6 shades and can be also used for either press-to-metal, or all-ceramic applications. They are usually pressed to full contour and then a dentin shade stain is applied and fired.
3. Neutral ingots are monochromatic and available in 6 types of ingots with different opacity/translucency. They are used for all-ceramic applications, and pressed to full contour and a dentin shade stain is applied and fired.

INSTRUCTIONS

PRESSED LAMINATE VENEERS (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

Ceramco iC ingots may be pressed as thin as 0.3 mm. Refer to the Ceramco iC All-Ceramic Directions for Use for procedures.

Section #6, Image #1: Waxed to Full Contour

Section #6, Image #2: Sprued

Section #6, Image #3: Pressed from a Translucent Ingot

Section #6, Image #4: Fitted to the Model

Section #6, Image #5: Mix the Dentin Shade Stain to a Creamy Consistency

Section #6, Image #6: Mix the Enamel Stain or Overglaze to a Creamy Consistency

Section #6, Image #7: Apply the Dentin Shade Stain

Section #6, Image #8: Vibrate to Smooth

Section #6, Image #9: Apply the Enamel Stain or Overglaze and Vibrate to Smooth

Section #6, Image #10: Completed Veneers

Section #6, Image #11: Facial Measurement

Section #6, Image #12: 0.3 mm Final Thickness

*Vita Lumin[®] and 3D-Master[®] are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co.

FELDSPATHIC LAMINATE VENEERS (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

Ceramco iC veneering porcelains may also be used to fabricate feldspathic veneers. They may be fabricated either on a platinum/palladium foil, or on a phosphate bonded investment that has a compatible CTE with the Ceramco iC porcelains.

MODEL PREPARATION – INVESTMENT TECHNIQUE

1. Use a super hard die stone to fabricate a model. White, ivory or buff colored die stones are preferred.
2. If it is necessary to mask out the underlying tooth structure, apply two coats of die spacer to within 1.0 mm of the prepared margin or to the base of the axial wall. This die spacer will create a 0.1 mm space for the resin cement. Also be sure to block out any undercuts.
3. Duplicate the gypsum model using Deguform® silicone duplicating material.
4. Mix the DENTSPLY Prosthetics Pressed Ceramic Investment Powder according to the Directions for Use. For each 100 gm packet of investment use the following ratio:
 - a. Crowns/veneers: 25 ml liquid
5. Allow to set a minimum of 30 minutes prior to removing from the Deguform mould. If necessary, construct a model with removable dies. Keep the models as small as possible to avoid long rehydration and dry/preheat cycles.

Section #7, Image #1: Refractory degassed

BURNOUT/DEGASSING THE INVESTMENT

Note: This procedure MUST be done prior to porcelain application. Failure to do so will result in porcelain porosities. We recommend using a burnout oven for this procedure. Gasses are given off when heating which could harm and contaminate a porcelain furnace.

Place the investment dies/model in a burnout oven heated to 850°C and hold for 30 min. Large models may require an additional hold time. At this point the investment should be completely white. Remove and allow to cool to room temperature.

PORCELAIN APPLICATION

Note: The first layer should never be opaque. Opaque does not completely etch and should never be used.

1. First soak the model and/or dies in distilled water for 5-10 minutes, or until no bubbles are coming to the surface.

Section #7, Image #2: Hydrate the refractory

2. If masking is desired, apply a thin coat of the Ceramco iC Opaque Dentin to the die. If no masking is desired, use a dentin or enamel porcelain and apply a thin layer to the die. Tap the die lightly and blot away the excess moisture. Fire the case according to the recommended firing cycle. This first layer will help to seal the die and aid in the application of the porcelain.

Section #7, Image #3: Apply a wash coat of porcelain

Section #7, Image #4: Wash coat fired

3. Resoak the model and/or die in distilled water for 5 - 10 minutes, or until no bubbles are coming to the surface.
4. Apply dentin and enamel porcelains to complete the case. Refer to the Ceramco iC PFM Directions for Use for specifics. Be sure to use the firing chart for refractory laminate veneers.

Section #7, Image #5: Apply the dentin and enamel

5. Finish, stain and glaze in your normal manner.
6. Once completed, remove the investment with aluminum oxide at no more than 2 bars (29 psi).

Section #7, Image #6: Completed veneers

Section #7, Image #7: Measuring the thickness

Section #7, Image #8: Measured thickness of 0.3 mm

7. If you need to do a repair after divesting, Ceramco iC Final Margin may be used. This material fires at a low temperature and will not distort your restoration.

ETCHING THE COMPLETED RESTORATION

The completed restoration must be etched prior to placement. This step is typically performed in the laboratory.

1. Use a steam cleaner to clean the inside of the completed crown or place the restorations in distilled water and clean in an ultrasonic cleaner for 10 minutes.
2. Apply a thin layer of DENTSPLY Prosthetics Etchant Gel to the surfaces of the restoration to be bonded to the prepared tooth structure. Apply with a disposable brush or plastic instrument.
3. Allow the etchant to stay in place for 1-2 min.
4. Remove the bulk of the etchant with a brush or cotton tip applicator.
5. Prepare the DENTSPLY Prosthetics neutralizer by mixing the package contents with 1 pint or 0.5 liter of water. Place the crown in the solution for approximately 20 seconds.

Note: Retain the remaining neutralizer mix in a plastic container with a secure lid. When the solution color changes from blue to yellow-green, you should discard it.

6. Clean the prepared surfaces with a steam cleaner or use distilled water in and ultrasonic cleaner for 10 minutes.
7. Do not touch the prepared crowns or place them on a model. This may contaminate the etched surface and decrease the bond strength.

Press-to-Metal Directions for Use

INDICATIONS FOR USE

The Ceramco iC Porcelain system is indicated for:

1. Porcelain Fused-to-Metal Applications: Single and multiple unit porcelain fused-to-metal fixed prosthodontic restorations.
2. Press-to-Metal Applications: Single and multiple unit press-to-metal fixed prosthodontic restorations.
3. All-Ceramic Applications: Single unit anterior and posterior premolar metal-free fixed prosthodontic restorations, laminate veneers, inlays and onlays.

CONTRAINDICATIONS

Only the indications listed above are suitable.

WARNINGS

Consult applicable DENTSPLY SDS (Safety Data Sheets) for safe handling and usage. Use only with recommended Ceramco iC system components. Use protective equipment to guard against dust when grinding.

PRECAUTIONS

Refer to the special notes throughout the instructions for best practices for using Ceramco iC materials. Other porcelains are not thermally compatible and will not have the correct dilatometric characteristics, and therefore should not be used. While initial results with some materials may appear acceptable, internal stress can compromise long term success.

ADVERSE REACTIONS

The Ceramco iC system contains components which may cause skin dryness, sensitization (allergic contact dermatitis or other allergic reactions in susceptible persons).

STORAGE

Store materials in packages with the lid tightly sealed. Do not leave open bottles where dust and debris may contaminate the porcelain. Replace cap when not in use. Avoid vibration and heat. Keep dry. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence. Do not use after expiration date.

INGOT TYPES

The Ceramco iC system offers three different types of shaded ingots.

1. Value Series ingots have a higher translucency, come in 6 shades and can be also used for either press-to-metal, or all-ceramic applications. They are usually pressed to full contour and then a dentin shade stain is applied and fired.

Note: For optimal efficiency and shade performance the Value Series and translucent ingots are recommended for press-to-metal restorations.

2. Neutral ingots are monochromatic and available in 6 types of ingots with different opacity/translucency. They are used for all-ceramic applications, and pressed to full contour and a dentin shade stain is applied and fired.
3. Dentin Shaded ingots come in all Vita Lumin[®], 3D-Master^{®*}, and bleached shades. They are slightly higher in opacity and generally used with an enamel add-on technique for either a press-to-metal, or all-ceramic application.

ALLOY SELECTION

Ceramco iC Press ingots are compatible with high noble, noble and predominately base alloys that are copper free and less than 10% silver. Be sure to consult the alloy manufacturer for alloy compositions and coefficient of thermal expansion data. Base Alloys should have a CTE of 13.8 – 14.1 @ 500°C. High noble and noble alloys should have a CTE of 14.0 - 14.5 @ 500°C.

INSTRUCTIONS

SUBSTRUCTURE DESIGN (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

1. Always follow the recommended substructure design for conventional porcelain fused-to-metal restorations.
 - a. Maintain the proper connector size.
 - b. Do not allow for more than 2.0 mm of unsupported porcelain.

Section #10, Image #1: Metal Coping

2. Porcelain margins may be pressed. However, these ingots are more translucent than other dentin porcelains. Unless the margins are prepared with a perfect shoulder or deep chamfer, the margins may

* Vita Lumin[®] and 3D-Master[®] are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co.

have too much translucency. A special Ceramco iC Margin porcelain may be utilized to assist in the masking of these margins. Refer to that section of these Directions for Use for specifics.

3. If the pontic is large, hollow pontics may be required. They will allow for a more rapid dissipation of the heat and help to prevent cracking.

PASTE OPAQUE APPLICATION

1. Oxidize or degass and clean the alloy according to the alloy manufacturer's instructions.
2. Check all castings to make sure that all surfaces are dry.
3. Use a flat brush and using a light touch apply an even coat of the paste opaque onto the metal substructure. Use Ceramco iC Opaque Modifier Liquid (not water) to adjust the paste consistency. Achieve an even coat and do not allow the material to puddle. Vibrate the restoration lightly to smooth the surface. If using a hollow pontic, fill in the hollow area to allow for a more even shade.

Note: Ceramco iC Opaque Crystals are NOT used for the press-to-metal technique.

4. Dry and fire according to the recommended temperatures. The correct visual indicator is an eggshell finish.
5. If necessary apply the second opaque porcelain layer using the same technique as for the first layer. Make sure this layer is thick enough to completely mask the metal substructure.
6. Dry and fire according to the recommended temperatures. Correct visual indicator is an eggshell finish.

Section #9, Image #1: Opaque with NO CRYSTALS

Section #10, Image #2: Opaque with NO CRYSTALS

CERAMCO IC PTM MARGIN PROCEDURES (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

The Ceramco iC Press Margin Porcelain is designed to assist in the esthetics of pressed porcelain margins by decreasing the translucent or gray margins. It is a single porcelain firing procedure followed by the normal waxing and pressing procedures.

METAL SUBSTRUCTURE

1. The crown should be prepared with a butt margin or deep chamfer. Design the metal substructure so the metal ends approximately 0.5 - 1.0 mm from the margin.

Section #8, Image #1: Coping with the prepared margin cut-back

2. Prepare and oxidize the alloy as recommended by the alloy manufacturer's directions for use.

OPAQUE APPLICATION

1. Apply and fire the Ceramco iC Paste using your normal procedures. Be sure the opaque completely covers the metal at the gingival margin. Failure to do so will make the completed porcelain margin appear gray.

Section #8, Image #2: Opaque with NO CRYSTALS

DIE SEALER APPLICATION

1. Using the DENTSPLY Prosthetics Die Sealer apply a thin, even layer of sealer over the die margin area.
2. Blow off any excess die sealer and allow the liquid to thoroughly dry.
3. Reapply the die sealer at least 2 more times and allow to dry.

DIE RELEASE APPLICATION

1. Using the DENTSPLY Prosthetics Die Release apply a thin, even layer of die release over the margin area. Blow off the excess die release.
2. Apply and allow to dry several more times until the die has a very slight sheen.

PORCELAIN MARGIN APPLICATION

1. Mix the margin porcelain with the Ceramco iC Margin and Opaque Corrector Liquid to a thick consistency. This mixture will air harden in approximately 5 minutes. If additional working time is required, mix in a small amount of distilled water.
2. Place the porcelain mixture in the gingival area, pushing the mixture down to but not over the gingival margin. Do not overbuild.
3. Allow the margin porcelain to dry completely (at least 5 minutes). The porcelain surface will be chalky and hard. Carefully lift the coping from the die.
4. Dry and fire the restoration according to the recommended temperatures. The correct visual indicator will be shiny yet grainy.

Section #8, Image #3: Margin porcelain applied and fired. Do not do a second application.

CROWN COMPLETION

Wax directly over the fired porcelain margin completing the proper contours. In most cases, the final margin will be made from the pressing process. This will increase the accuracy of the porcelain margin, and eliminates the second margin firing. If after staining and glazing, a low-temperature margin add-on is required, use the Ceramco iC Final Margin. Refer to the Ceramco iC PFM Directions for Use for specifics.

Section #8, Image #4 Full contour wax-up

Section #8, Image #5: Wax over the margin

Section #8, Image #6: Completed crown

DETERMINING WHICH INGOTS TO USE

The Ceramco iC system has two different shading formulations. Dentin Shaded ingots and Value Series ingots. The 46 dentin shaded ingots are designed to be pressed and then veneered with Ceramco iC Enamel porcelains. There are six Value Series ingots that are designed to be pressed to full contour, and then stained to achieve the final shade.

DETERMINING HOW MANY INGOTS TO USE

Ingot Size	Ingot Diameter	Ingot Type
2 gram	10mm	Dentin Shaded, Value Series
2 gram	15mm	Value Series
5 gram	15mm	Value Series

The following combination of rings and ingots may be used:

Ring Size	# of Ingots	Ingot Sizes
100 g	1 or 2	2 gram, 10 mm diameter
200 g	1 or 2	2 gram, 10 mm diameter
300 g	1 or 2	5 gram, 15 mm diameter
	1 – 5	2 gram, 15 mm diameter
	1 Up to 2	5 gram, 15 mm diameter 2 gram, 15 mm diameter

1. Weigh the sprue base and the opaqued metal substructure.
2. Wax the case and then attach it to the sprue base and reweigh the sprue base.
3. The difference between the two is the weight of the wax. For 4-unit bridges and higher add 2g.
4. Use the following table for ingot size and selection.

Wax Weight	Ingot Size	# Ingots
0.0 – 0.60 grams	2 gram	1
0.61 – 1.40 grams	2 gram 5 gram	2 1
1.41 – 2.5 grams	5 gram	2

FULL CONTOUR WAX-UP PROCEDURES

1. Wax the crown using a conventional inlay wax. Note: Use a non-contaminating inlay wax. After burning out, some inlay waxes may leave an ash or carbon residue. This will result in dark areas on the pressed crown.
2. Recreate all necessary anatomical features in wax.
3. To achieve the optimum shade, the recommended wax thickness is 1.2 mm. The minimum thickness of the wax is 0.8 mm.

Section #10, Image #3: Wax to Full-Contour

Section #9, Image #2: Wax to Full-Contour

SPRUING PROCEDURES

1. All sprues should be 2-3 mm in diameter or 8 gauge/3mm round. Extremely large crowns may require a 6 gauge/4mm round sprue.
2. Sprues should be approximately 2-4 mm in length.
3. For anterior teeth, attach a single sprue directly to each incisal edge. For posterior teeth, attach a single sprue at a 45° on the largest cusp tip. Do not use sharp angles.

Section #10, Image #4: Sprued

4. Attach the sprues to the sprue base in the same manner. Slightly flare the wax at the point of attachment.
5. Allow a minimum of 3.0 mm between each wax-up, and 8.0 mm between the wax-ups and the exterior wall of the investment.
6. Allow a minimum of 8.0 mm between the top of the wax-up and the leveling ring.

INVESTING PROCEDURES

Note: You must use a phosphate bonded investment such as DENTSPLY Prosthetics PRESSED CERAMIC INVESTMENT POWDER and LIQUID for these procedures.

1. For each 100 gram packet of investment use 20 ml of investment liquid and 5 ml of distilled water.
2. Place the liquid in a slightly moistened mixing bowl followed by the investment.
3. Hand-mix the materials for 15 seconds.
4. Use a full vacuum and mix the materials under slow speed for 60 seconds.
5. Carefully paint the restorations with the investment and pour the remaining investment into the silicone ring.
6. Place the leveling ring on the top of the silicone ring. Excess investment will push through the top of the leveling ring. Be sure it is fully seated. This step is important, as it will determine the proper height and angle of the investment for the pressing furnace.

BURNOUT PROCEDURES

1. Allow the investment to bench set a minimum of 20 minutes.
2. Remove the leveling ring and the sprue base from the investment.
3. For the small 100 gram and large 200 gram rings:
 - a. Use a knife to carefully remove any rough investment from the top of the ring.

- b. The top and bottom of the ring should be flat. Do not use a model trimmer. Do not decrease the height of the ring.
- c. Be sure to maintain the 90° angle between the ends and the sides of the ring.

- 4. For the extra large 300 gram ring:
 - a. Use a knife to carefully remove the rough investment from the leveling ring.
 - b. DO NOT REMOVE THE RAISED INVESTMENT ON THE BASE. This is used to position the ring in the pressing furnace.
- 5. Do not allow any investment or debris to enter the sprue hole. If the ring is allowed to set over one hour, place the ring in a humidior until ready to burnout.

Note: Be sure the floor of the burnout furnace is clean. Small pieces of debris will enter the sprue hole and will result in small black specks in the pressed ceramic. Be sure the floor of the burnout furnace has ridges to allow the wax to easily escape. Placing the ring in the burnout furnace on a flat surface may cause the ring to crack.

OPTION #1: THIS IS THE PREFERRED TECHNIQUE. USED WHEN THE INVESTMENT HAS BENCH SET A MINIMUM OF 20 MINUTES

- 1. If desired, preheat the burnout furnace to 850°C/1562° F.
- 2. Do not pre-soak the investment in water prior to placing in the burnout furnace.
- 3. Place the rings (sprue hole down) directly into the preheated furnace.
- 4. Place the alumina plungers into the preheated furnace.

Note: If disposable refractory plungers are used do not preheat.

- 5. Heat soak the 100g rings 45 minutes, the 200g rings 60 minutes, and the 300g rings 90 minutes.
- 6. Add 15 minutes for each additional ring.

OPTION #2: USED WHEN THE INVESTMENT HAS BENCH SET OVER 60 MINUTES

- 1. If the investment has become dry, rehydrate the investment by placing the ring in a small amount of water for 10-15 minutes prior to burnout.
- 2. Place the rings and plungers into a burnout furnace at room temperature.

Note: If disposable refractory plungers are used do not preheat.

- 3. Raise the temperature of the furnace at 15°C/27°F per minute to 850°C/1562°F.

- 4. Heat soak the 100g rings 45 minutes, the 200g rings 60 minutes, and the 300g rings 90 minutes.
- 5. Add 15 minutes for each additional ring.

PRESSING PROCEDURES

- 1. The program is as follows:

Low Temperature:	700°C/1292°F
Heat Rate:	60°C/108°F per min
High Temperature:	890°C/1634°F
Hold Time:	20 min
Pressing Time:	20 min (100g ring)
	20 min (200g ring)
Cool Time:	0 min
- 2. Remove a ring from the burnout furnace and place it on a heat resistant surface with the sprue hole up.
- 3. Carefully place the correct ingot(s) into the sprue hole.
- 4. Remove one of the preheated alumina plungers from the burnout furnace or use a disposable refractory plunger and place it in the sprue hole on top of the ingots.
- 5. Place the ring in the center of the firing platform of the pressing furnace and press the Start button. Use pressing chart on page 25.
- 6. After the pressing cycle has been completed, remove the ring from the furnace and allow to bench set to cool. The ring may be fast cooled by placing it in front of a circulating fan.

DIVESTING PROCEDURES

- 1. Take an alumina plunger or unused refractory plunger and place alongside the pressed ring.
- 2. Use a pencil and mark the approximate length of the used plunger.
- 3. If using an alumina plunger use a disk or 50-micron aluminum oxide and cut the investment down to the top of the plunger. Carefully remove the plunger from the investment. If the plunger is thoroughly cleaned with aluminum oxide it may be used again. If any debris remains on the plunger and it is used again, it will harm the investment and will result in a mispress.
- 4. Use 50-micron aluminum oxide with 60 psi (4.14 bars) to remove the bulk of the investment. Do not blast the pressed ceramic at this pressure with the aluminum oxide. The aluminum oxide is too abrasive and will harm the pressed ceramic.
- 5. Use 50-micron glass beads with 20 psi (1.38 bars) to remove the investment around the pressed ceramic. Be sure to remove all the investment.

Section #10, Image #5: Pressed in a Dentin Shaded ingot

Section #9, Image #3: Pressed in a Value Series Ingot

FINISHING PROCEDURES

1. Use a thin diamond disk to remove the sprues from the ceramic crowns.
2. Use a diamond bur or suitable aluminum oxide stones to recontour the sprue attachment. Do not use abrasives that create excess heat, such as worn carbides or coarse grit stones.
3. Correct contours as necessary.

VALUE SERIES INGOTS STAINING AND GLAZING (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

1. Clean the porcelain surface by lightly blasting with 50-micron aluminum oxide at 20 psi, followed by cleaning with a steam cleaner or in an ultrasonic cleaner with distilled water.
2. Place a small amount of the Dentin Shade Stain on the pallet. Mix in the Stain and Glaze Liquid to achieve a creamy viscosity. Apply a liberal amount of the mixture to the porcelain surface. Vibrate lightly to smooth out the shade stain. Apply the stain to the dentin area until the desired effect is achieved.

Section #9, Image #4: Mix the Dentin Shade Stain to a Creamy Consistency

Section #9, Image #6: Apply the Dentin Shade Stain

Section #9, Image #7: Vibrate to smooth

3. Place a small amount of the Enamel Stain or Overglaze on the pallet. Mix in the Stain and Glaze Liquid to achieve a creamy viscosity and apply the mixture only to the incisal edge surface. Feather along the dentin-incisal interface. Vibrate lightly to smooth out the stain/overglaze.

Section #9, Image #5: Mix the Enamel Stain or Overglaze to a Creamy Consistency

Section #9, Image #8: Apply Enamel Stain or Overglaze to the Surface

4. The Ceramco iC colored stains may be used for individual characterizations.

Section #9, Image #9: Apply Stains to the Occlusal Surface

5. Fire according to recommended temperatures.

Section #9, Image #10: Completed Crown

Section #9, Image #11: Completed Crown

DENTIN SHADED INGOTS ENAMEL/ADD-ON PORCELAIN (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

Any of the Ceramco iC veneering porcelains may be used for any add-on requirements. Follow the Ceramco iC Porcelain Fused-to-Metal Directions for Use for any add-on procedures. Be sure the porcelain surfaces are clean prior to any porcelain application.

Section #10, Image #6: Enamel cut-back

Section #10, Image #7: Enamel application

Section #10, Image #8: Completed Crown

All-Ceramic Directions for Use

INDICATIONS FOR USE

The Ceramco iC Porcelain system is indicated for:

1. Porcelain Fused-to-Metal Applications: Single and multiple unit porcelain fused-to-metal fixed prosthodontic restorations.
2. Press-to-Metal Applications: Single and multiple unit press-to-metal fixed prosthodontic restorations.
3. All-Ceramic Applications: Single unit anterior and posterior premolar metal-free fixed prosthodontic restorations, laminate veneers, inlays and onlays.

CONTRAINDICATIONS

Only the indications listed above are suitable.

WARNINGS

Consult applicable DENTSPLY SDS Safety Data Sheets) for safe handling and usage. Use only with recommended Ceramco iC system components. Use protective equipment to guard against dust when grinding.

PRECAUTIONS

Refer to the special notes throughout the instructions for best practices for using Ceramco iC materials. Other porcelains are not thermally compatible and will not have the correct dilatometric characteristics, and therefore should not be used. While initial results with some materials may appear acceptable, internal stress can compromise long term success.

ADVERSE REACTIONS

The Ceramco iC system contains components which may cause skin dryness, sensitization (allergic contact dermatitis or other allergic reactions in susceptible persons).

STORAGE

Store materials in packages with the lid tightly sealed. Do not leave open bottles where dust and debris may contaminate the porcelain. Replace cap when not in use. Avoid vibration and heat. Keep dry. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence. Do not use after expiration date.

INGOT TYPES

The Ceramco iC system offers three different types of shaded ingots.

1. Dentin Shaded ingots come in all Vita Lumin[®], 3D-Master^{®*}, and bleached shades. They are slightly higher in opacity and generally used with an enamel add-on technique for either a press-to-metal, or all-ceramic applications.
2. Value Series ingots have a higher translucency, come in 6 shades and can be also used for either press-to-metal, or all-ceramic applications. They are usually pressed to full contour and then a dentin shade stain is applied and fired.
3. Neutral ingots are monochromatic and available in 6 types of ingots with different opacity/translucency. They are used for all-ceramic applications, and pressed to full contour and a dentin shade stain is applied and fired.

INSTRUCTIONS

MODEL PREPARATION

1. Use a super hard die stone to fabricate a cast with removable dies. White, ivory or buff colored die stones are preferred. These die stones more closely replicate the color of natural dentition and will make it easier to determine the accuracy of the final shade.
2. If it is necessary to mask out the underlying tooth structure, apply two coats of die spacer to within 1.0 mm of the prepared margin or to the base of the axial wall. This die spacer will create a 0.1 mm space for the resin cement.
3. Allow the die spacer to completely harden.

WAX-UP CONSIDERATIONS

ENAMEL LAYERED TECHNIQUE:

This technique offers excellent esthetics and high productivity. Wax a crown to full contour, and cutback only the enamel area. This crown is pressed using one of the more opaceous dentin shaded ingots. The Ceramco iC Enamel porcelains are then used to complete the build-up.

STAINING TECHNIQUE:

This technique is highly productive and due to the thicker core material the result is a more translucent

*Vita Lumin[®] and 3D-Master[®] are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co.

restoration with a higher strength. This technique is recommended for premolars and anterior restorations that are subjected to high occlusal stress. Wax the crown to full contour. This crown is pressed using one of the Ceramco iC Value Series (VS) or Neutral (N) ingots and completed using the Ceramco iC Dentin Shade Stain porcelains.

WAX-UP PROCEDURES (REFER TO THE INSIDE COVER FOR VISUAL INSTRUCTION)

1. Apply a thin coat of a conventional die lubricant.
2. Wax the crown using a conventional inlay wax.

Note: Use a non-contaminating inlay wax. After burning out, some inlay waxes may leave an ash or carbon residue. This will result in dark areas on the pressed crown.

3. Recreate all necessary anatomical features in wax and completely seal the margins. For full coverage restorations, be sure the wax has a minimum thickness of 0.8 mm. Laminate veneers should have a minimum thickness of 0.3 mm.

Section #11, Image #1: Wax to full contour

Section #12, Image #1: Wax to full contour

DETERMINING HOW MANY INGOTS TO USE

Ingot Size	Ingot Diameter	Ingot Type
2 gram	10mm	Dentin Shaded, Value Series
2 gram	15mm	Value Series
5 gram	15mm	Value Series

The following combination of rings and ingots may be used:

Ring Size	# of Ingots	Ingot Sizes
100 g	1 or 2	2 gram, 10 mm diameter
200 g	1 or 2	2 gram, 10 mm diameter
300 g	1 or 2	5 gram, 15 mm diameter
	1 – 5	2 gram, 15 mm diameter
	1 Up to 2	5 gram, 15 mm diameter 2 gram, 15 mm diameter

1. Weigh the sprue base and the opaqued metal substructure.
2. Wax the case and then attach it to the sprue base and reweigh the sprue base.
3. The difference between the two is the weight of the wax. For 4-unit bridges and higher add 2g.
4. Use the following table for ingot size and selection.

Wax Weight	Ingot Size	# Ingots
0.0 – 0.60 grams	2 gram	1
0.61 – 1.40 grams	2 gram 5 gram	2 1
1.41 – 2.5 grams	5 gram	2

ATTACHING THE SPRUES

1. All sprues should be 2-3 mm in diameter or 8 gauge/3 mm round. Extremely large crowns may require a 6 gauge/4 mm round sprue.
2. Sprues should be 2-4 mm in length.
3. Attach a single sprue directly to the incisal edge or on a cusp tip. Placing the sprue at an angle to the incisal edge may result in the investment fracturing during the pressing procedures.
4. If a thin area of wax separates two thick areas of wax, a sprue should be attached to each of the thick areas.
5. Do not flare the wax at the point of attachment. Do not use sharp angles.
6. After removing the wax-up from the die, be sure the die spacer did not adhere to the inside of the crown.
7. Attach the sprues to the sprue base in the same manner. Slightly flare the wax at the point of attachment. Do not use sharp angles.

Section #11, Image #2: Sprued

8. Allow a minimum of 3.0 mm between each wax-up, and between the wax-ups and the exterior wall of the investment.
9. Allow a minimum of 8.0 mm between the top of the wax-up and the leveling ring.

INVESTING PROCEDURES

Note: You must use a phosphate bonded investment such as DENTSPLY Prosthetics Pressed Ceramic Investment Powder and Liquid for these procedures.

1. Debubblers may be used, but ensure they have completely dried on the wax patterns. If the wax patterns are still wet, the debubbler may leave a residue that will result in a rough surface on the pressed ceramic. Wax pattern cleaners may be used, but no residue should remain.
2. For each 100 gm packet of investment use the following ratios:
 - a. Crowns/veneers: 20ml liquid 5ml distilled water
 - b. Inlays/onlays: 12.5ml liquid 12.5ml distilled water

Note: The suggested liquid/water ratio will normally give you the correct expansion. The investment liquid controls the expansion of the investment. Using more investment liquid and less distilled water will result in more expansion. Using less investment liquid and more distilled water will result in less expansion. Do not change the powder/liquid ratio, as this will affect the strength and the working time of the investment. Do not dilute the liquid/distilled water ratio below 50%.

3. Place the liquid in a slightly moistened mixing bowl followed by the investment.
4. Hand-mix the materials for 15 seconds.
5. Use a full vacuum and mix the materials under slow speed for 60 seconds.
6. Carefully paint the restorations with the investment and pour the remaining investment into the ring.

Section #11, Image #3: Invest

7. Place the leveling ring on top of the silicone ring. Excess investment will push through the top of the leveling ring. Be sure it is fully seated on the silicone ring. This step is important, as it will determine the proper height and angle of the investment for the pressing furnace.

BURNOUT PROCEDURES

PREPARING THE INVESTED RING FOR BURNOUT

1. Allow the investment to bench set a minimum of 20 minutes.
2. Remove the silicone ring, leveling ring and the sprue base from the investment.
3. Using a knife, carefully remove any rough investment from the top and bottom of the ring. The top and bottom of the ring should be flat. Do not use a model trimmer. Do not decrease the height of the ring. Be sure to maintain the 90° angle between the ends and the sides of the ring.
4. Do not allow any investment or debris to enter the sprue hole. If the ring is allowed to set over one hour, place the ring in a humidior until ready to burnout.

Note: Be sure the floor of the burnout furnace is clean. Small pieces of debris will enter the sprue hole and will result in small black specks in the pressed ceramic. Be sure the floor of the burnout furnace has ridges to allow the wax to easily escape. Placing the ring in the burnout furnace on a flat surface may cause the ring to crack. If possible, utilize a separate burnout furnace only for all-ceramic restorations.

OPTION #1 (PREFERRED): USED WHEN THE INVESTMENT HAS BENCH SET ONLY 20 MINUTES

1. If desired, preheat the burnout furnace to 850°C/1562° F or raise the temperature of the furnace at 15°C/27°F per minute to 850°C/ 1562°F.
2. Do not pre-soak the investment in water prior to placing in the burnout furnace.
3. Place the rings (sprue hole down) directly into the preheated furnace.
4. Place the alumina plungers into the preheated furnace.

Note: If disposable refractory plungers are used do not preheat.

5. Allow the small rings to heat soak 45 minutes and the large rings 60 minutes.
6. Add 15 minutes for each additional ring. Additional time may be required for large rings.

OPTION #2: USED WHEN THE INVESTMENT HAS BENCH SET OVER 60 MINUTES

1. If the investment has become dry, rehydrate the investment by placing the ring in water for 10-15 minutes prior to burnout.
2. Place the rings and plungers into a burnout furnace at room temperature.

Note: If disposable refractory plungers are used do not preheat.

3. Be sure the floor of the burnout furnace has ridges to allow the wax to easily escape. Placing the ring in the burnout furnace on a flat surface may cause the ring to crack.
4. Raise the temperature of the furnace at 15°C/27°F per min to 850°C/1562°F.
5. Allow the small rings to heat soak 45 min and the large ones 60 min.
6. Add 15 minutes for each additional ring.

PRESSING PROCEDURES

1. The program is as follows:

Low Temperature:	700°C/1292°F
Heat Rate:	60°C/108°F per min
High Temperature:	890°C/1634°F
Hold Time:	20 min
Pressing Time:	20 min (100g ring) 20 min (200g ring)
Cool Time:	0 min
2. Remove a ring from the burnout furnace and place it on a heat resistant surface with the sprue hole up.
3. Carefully place the correct ingot(s) into the sprue hole.
4. Remove one of the preheated alumina plungers from the burnout furnace or use a disposable refractory plunger and place it in the sprue hole on top of the ingots.
5. Place the ring in the center of the firing platform of the pressing furnace and press the Start button. Use pressing chart on page 25.
6. After the pressing cycle has been completed, remove the ring from the furnace and allow to bench set to cool. The ring may be fast cooled by placing it in front of a circulating fan.

DIVESTING PROCEDURES

1. Take a plunger and place alongside the pressed ring.
2. Use a pencil and mark the approximate length of the used plunger.
3. If using an alumina plunger use a disk or 50-micron aluminum oxide and cut the investment down to the top of the plunger. Carefully remove the alumina plunger from the investment. If the plunger is thoroughly cleaned with aluminum oxide it may be used again. If any debris remains on the plunger and it is used again, it will harm the investment and will result in a mispress.
4. Use 50-micron aluminum oxide with 60 psi (4.14 bars) to remove the bulk of the investment. Do not blast the pressed ceramic at this pressure with the aluminum oxide. The aluminum oxide is too abrasive and will harm the pressed ceramic.
5. Use 50-micron glass beads with 20 psi (1.38 bars) to remove the investment around the pressed ceramic. Be sure to remove all the investment.

Section #11, Image #4: Press using a Dentin Shaded Ingot

Section #12, Image #2: Pressed

FINISHING PROCEDURES

1. Conventional abrasives used for dental porcelains may be used for the Ceramco iC All-Ceramic. Use a thin diamond disk to remove the sprues from the ceramic crowns.
2. Use a diamond bur or suitable aluminum oxide stones to recontour the sprue attachment. Do not use abrasives that create excess heat, such as worn carbides or course grit stones. This will cause a thermal crack and the ceramic coping/crown will have to be remade.

Section #11, Image #5: Enamel Cut-Back

Section #12, Image #3: Seated on the Model

3. Closely examine the inside of the restoration for bubbles or irregularities. These can be removed with a thin diamond bur or stone.
4. Carefully seat the restoration on the die. Indicating materials may be used to aid in seating the restoration. Be sure indicating materials have been completely removed prior to placing the restorations in the porcelain furnace.
5. Use a gray silicone wheel or point to adjust the margins. Other abrasives may chip the pressed ceramic and decrease the marginal integrity.

Note: *If the margins are short, Ceramco iC Dentin porcelains or the Ceramco iC Final Margin may be used to perfect the*

marginal integrity during the dentin or glaze firing. The regular Ceramco iC Margin porcelains are too opaque and not recommended for use.

COMPOSITE DIE PROCEDURES

The DENTSPLY Prosthetics Composite Die Material is designed to mimic the actual shade of the patients prepared tooth. When this material is placed inside the pressed ceramic core, it will assist in accurate shade reproductions.

If the dentist did not take the shade of the prepared tooth, the composite die material may be used to verify to final shade. Select the proper composite die material from the appropriate table.

1. Apply the DENTSPLY Prosthetics Die Release to the inside of the ceramic restoration and allow it to dry.
2. Place a small amount of the DENTSPLY Prosthetics Composite Die Material on the inside of the restoration. Pack the material to remove any voids. Immediately push a dowel pin into the uncured composite die material. Remove any excess composite from the margin area.
3. Light cure the composite for 1-2 minutes using a hand held light curing unit or the Triad® 2000 curing unit (DENTSPLY/Trubyte).
4. Remove the composite die material from the restoration and clean it using a steam cleaner or in distilled water in an ultrasonic cleaner for 10 minutes.

PORCELAIN APPLICATION

Note: *For Ceramco iC All-Ceramic restorations, fill inside of the restoration with heat resistant refractory material (SuperPeg™II* Refractory Material available from DENTSPLY Prosthetics), then insert DENTSPLY Prosthetics All-Ceramic firing pins, and place the restoration on a honeycomb sagger tray to fire. DENTSPLY Prosthetics pins are recommended since they are easy to bend and do not act as a heat sink. Other firing pins will not allow the proper heat transfer and will create a stress point. While initial results with some firing pins may appear acceptable, internal stress can compromise long term success.*

ENAMEL LAYERED TECHNIQUE

1. Use 50-micron aluminum oxide at 20 psi and lightly blast the exterior surface of the ceramic restoration. Be careful not to harm the margins.
2. Use a steam cleaner to clean the surfaces or put the restoration in distilled water and place in an ultrasonic cleaner for 10 minutes.
3. Apply and fire the enamel porcelains to complete the restorations. All the components of the Ceramco iC system may be used to enhance the esthetics and contours of the restoration.

Section #11, Image #6: Enamel Effects Porcelain

Section #11, Image #7: Clear Porcelain

Section #11, Image #8: Enamel Completed

4. Always use the honeycomb sagger trays with the DENTSPLY Prosthetics All-Ceramic firing pins. Do not use other types of firing pins.
5. Refer to the Ceramco iC instructions for specific details and firing instructions.
6. The final shade may be checked using the DENTSPLY Prosthetics Composite Die Material.
7. To glaze the restoration, mix the Overglaze with the Stain and Glaze Liquid to a thin consistency and apply to the porcelain surface. Note the glazing cycle has a 2 minute hold time.

Section #11, Image #9: Completed Crowns

STAINING TECHNIQUE

1. Use 50-micron aluminum oxide at 20 psi and lightly blast the exterior surface of the crown. Be careful not to harm the margins.
2. Use a steam cleaner to clean the surfaces or put the restoration in distilled water and place in an ultrasonic cleaner for 10 minutes.
3. Always use the honeycomb sagger trays with the DENTSPLY Prosthetics All-Ceramic firing pins.
4. If necessary, corrections may be made to the crown using Ceramco iC porcelains.
5. Place the pressed crown on the prepared DENTSPLY Prosthetics Composite Die Material.
6. Place a small amount of the Dentin Shade Stain on the pallet. Mix in the Stain and Glaze Liquid to achieve a creamy viscosity. Apply a liberal amount of the mixture to the porcelain surface. Vibrate lightly to smooth out the shade stain. Apply the stain to the dentin area only until the desired effect is achieved.

Section #12, Image #4: Mix the Dentin Shade Stain to a Creamy Consistency

Section #12, Image #6: Apply the Dentin Shade Stain

Section #12, Image #7: Vibrate to Smooth out the Stain

7. Place a small amount of the Enamel Stain or Overglaze on the pallet. Mix in the Stain and Glaze Liquid to achieve a creamy viscosity and apply the mixture only to the incisal edge surface. Feather along the dentin-incisal interface. Vibrate lightly to smooth out the stain/overglaze.

Section #12, Image #5: Mix the Enamel Stain or Overglaze to a Creamy Consistency

Section #12, Image #8: Apply the Enamel Stain

8. The Ceramco iC colored stains may be used for individual characterizations.
9. Remove the die material from the restoration. To ensure proper fit, remove excess staining porcelain from the interior of the crown as well as the margin areas.
10. Fire the crown according to recommended firing cycles.

Note: Overglaze/Shade Stain/Color Stain - nominal high temperature 750°C and hold time 1.5 min. Temperature and hold time can be increased up to 760°C and 2.5 min., respectively. For furnaces that underfire, temperature can be increased by additional 10°C. Vacuum is recommended from low temperature of 500°C to high temp. Hold shall be in air. Heat rate is 35°C/min.

Note: The glazing cycle shall have hold time from 0.5 min. to 2.5 min. Shorter times are recommended for two applications of shade stain/overglaze.

Note: If higher sheen is desirable, firing temperature can be increased 10°C.

Section #12, Image #9: Completed Crowns

ETCHING THE COMPLETED RESTORATION

The completed restoration must be etched prior to placement. This step is typically performed in the laboratory.

1. Use a steam cleaner to clean the inside of the completed crown or place the restorations in distilled water and clean in an ultrasonic cleaner for 10 minutes.
2. Apply a thin layer of Ceramco iC Etchant Gel to the surfaces of the restoration to be bonded to the prepared tooth structure. Apply with a disposable brush or plastic instrument.
3. Allow the etchant to stay in place for 1-2 min.
4. Remove the bulk of the etchant with a brush or cotton tip applicator.
5. Prepare the DENTSPLY Prosthetics neutralizer by mixing the package contents with 1 pint or 0.5 liter of water. Place the crown in the solution for approximately 20 seconds.

Note: Retain the remaining neutralizer mix in a plastic container with a secure lid. When the solution color changes from blue to yellow-green, you should discard it.

6. Clean the prepared surfaces with a steam cleaner or use distilled water in an ultrasonic cleaner for 10 minutes.
7. Do not touch the prepared crowns or place them on a model. This may contaminate the etched surface and decrease the bond strength.

CERAMCO IC PFM PORCELAIN AND CUT-BACK INGOT CHART

Porcelain Fused-to-Metal						Press-to-Metal & All-Ceramic Cut-Back			
Ceramco iC Shade	Keyed to Vita® 3D-Master® ¹ Shade ²	Ceramco iC Enamel Porcelain Fused-to-Metal	High Chroma Dentin	Porcelain Fused-to-Metal Margin	Final Margin	Composite Die	Cut-Back Ingot Shade	Enamel	Margins Press-to-Metal
A1		Extra Light	HC2	M4	FM3	F1	A1	Extra Light	PM3
A2		Extra Light	HC2	M6	FM4	F12	A2	Extra Light	PM3
A3		Light	HC2	M7	FM6	F10	A3	Light	PM4
A3.5		Medium	HC5	M8	FM6	F9	A3.5	Medium	PM4
A4		Medium	HC5	M11	FM8	F7	A4	Medium	PM5
B1		Extra Light	HC1	M3	FM3	F11	B1	Extra Light	PM6
B2		Light	HC1	M5	FM4	F11	B2	Light	PM6
B3		Light	HC4	M8	FM6	F8	B3	Light	PM4
B4		Light	HC4	M8	FM6	F8	B4	Light	PM5
C1		Light	HC3	M9	FM5	F3	C1	Light	PM7
C2		Medium	HC3	M9	FM5	F4	C2	Medium	PM7
C3		Medium	HC6	M10	FM7	F5	C3	Medium	PM8
C4		Medium	HC6	M11	FM7	F6	C4	Medium	PM8
D2		Light	HC2	M8	FM6	F2	D2	Light	PM7
D3		Medium	HC1	M7	FM6	F2	D3	Medium	PM4
D4		Medium	HC6	M10	FM4	F3	D4	Medium	PM5
CC-000*	i7**	White	HC1	M1	FM1	F1	CC-000	White	PM2
CC-001	0M1	White	HC1	M1	FM1	F1	CC-001	White	PM2
CC-002	0M2	White	HC1	M2	FM2	F1	CC-002	White	PM1
CC-003	0M3	White	HC1	M2	FM2	F1	CC-003	White	PM1
CC-101	1M1	Extra Light	HC1	M2	FM2	F1	CC-101	Extra Light	PM1
CC-102	1M2	Extra Light	HC2	M4	FM3	F1	CC-102	Extra Light	PM1
CC-203	2L1.5	Extra Light	HC1	M5	FM4	F11	CC-203	Extra Light	PM6
CC-204	2L2.5	Extra Light	HC2	M7	FM6	F10	CC-204	Extra Light	PM6
CC-205	2M1	Extra Light	HC3	M6	FM5	F3	CC-310	Extra Light	PM6
CC-206	2M2	Extra Light	HC2	M4	FM3	F12	CC-206	Extra Light	PM6
CC-207	2M3	Extra Light	HC2	M7	FM6	F10	CC-207	Extra Light	PM6
CC-208	2R1.5	Extra Light	HC2	M6	FM5	F2	CC-208	Extra Light	PM6
CC-209	2R2.5	Extra Light	HC2	M4	FM3	F10	CC-209	Extra Light	PM6
CC-310	3L1.5	Extra Light	HC6	M10	FM7	F4	CC-310	Extra Light	PM7
CC-311	3L2.5	Light	HC4	M8	FM6	F9	CC-311	Light	PM8
CC-312	3M1	Extra Light	HC3	M6	FM5	F3	CC-312	Extra Light	PM6
CC-313	3M2	Light	HC6	M7	FM6	F4	CC-313	Light	PM4
CC-314	3M3	Light	HC5	M9	FM6	F9	CC-314	Light	PM8
CC-315	3R1.5	Light	HC4	M8	FM6	F3	CC-315	Light	PM7
CC-316	3R2.5	Light	HC5	M9	FM6	F7	CC-316	Light	PM8
CC-417	4L1.5	Light	HC6	M10	FM7	F4	CC-417	Light	PM4
CC-418	4L2.5	Light	HC5	M9	FM8	F7	CC-418	Light	PM5
CC-419	4M1	Light	HC6	M7	FM6	F3	CC-419	Light	PM7
CC-420	4M2	Light	HC5	M9	FM8	F9	CC-420	Light	PM8
CC-421	4M3	Light	HC5	M11	FM8	F8	CC-421	Light	PM8
CC-422	4R1.5	Light	HC5	M9	FM6	F5	CC-422	Light	PM7
CC-423	4R2.5	Light	HC5	M11	FM8	F9	CC-423	Light	PM8
CC-524	5M1	Medium	HC6	M10	FM7	F6	CC-524	Medium	PM8
CC-525	5M2	Medium	HC5	M11	FM8	F7	CC-525	Medium	PM5
CC-526	5M3	Medium	HC5	M11	FM8	F7	CC-526	Medium	PM8

* CC-000 is keyed to DENTSPLY illumine bleached shade i7.

** i7 is not a Vita Shade.

¹Vita® and 3D-Master® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. ²Data on file.

CERAMCO IC INGOT CHART - FULL CONTOUR SHADING

Shade	Keyed to Vita® 3D-Master®¹ Shade²	Composite Die	Ceramco iC Value Series Ingots	Dentin Shade Stains For Value Series and Neutral Ingots	Overglaze & Enamel Shade Stain[±]
A1		F1	VS1	A1	Overglaze
A2		F12	VS2	A2	Overglaze
A3		F10	VS2	A3	Overglaze
A3.5		F9	VS5	A3.5	Medium
A4		F7	VS5	A4	Medium
B1		F11	VS1	B1	Overglaze
B2		F11	VS1	B2	Overglaze
B3		F8	VS4	B3	Extra Light
B4		F8	VS4	B4	Light
C1		F3	VS3	C1	Overglaze
C2		F4	VS3	C2	Overglaze
C3		F5	VS3	C3	Overglaze
C4		F6	VS3	C4	Overglaze
D2		F2	VS3	D2	Overglaze
D3		F2	VS2	D3	Overglaze
D4		F3	VS3	D4	Overglaze
CC-000*	i7**	F1	VS0	CC-000	Overglaze
CC-001	0M1	F1	VS0	CC-001	Overglaze
CC-002	0M2	F1	VS0	CC-002	Overglaze
CC-003	0M3	F1	VS0	CC-003	Overglaze
CC-101	1M1	F1	VS1	CC-101	Overglaze
CC-102	1M2	F1	VS1	CC-102	Overglaze
CC-203	2L1.5	F11	VS1	CC-203	Overglaze
CC-204	2L2.5	F10	VS1	CC-204	Overglaze
CC-205	2M1	F3	VS1	CC-205	Overglaze
CC-206	2M2	F12	VS1	CC-206	Overglaze
CC-207	2M3	F10	VS4	CC-207	Extra Light
CC-208	2R1.5	F2	VS2	CC-208	Overglaze
CC-209	2R2.5	F10	VS2	CC-209	Overglaze
CC-310	3L1.5	F4	VS2	CC-310	Overglaze
CC-311	3L2.5	F9	VS4	CC-311	Light
CC-312	3M1	F3	VS3	CC-312	Overglaze
CC-313	3M2	F4	VS2	CC-313	Overglaze
CC-314	3M3	F9	VS4	CC-314	Light
CC-315	3R1.5	F3	VS2	CC-315	Overglaze
CC-316	3R2.5	F7	VS4	CC-316	Light
CC-417	4L1.5	F4	VS3	CC-417	Overglaze
CC-418	4L2.5	F7	VS5	CC-418	Medium
CC-419	4M1	F3	VS3	CC-419	Overglaze
CC-420	4M2	F9	VS5	CC-420	Light
CC-421	4M3	F8	VS5	CC-421	Medium
CC-422	4R1.5	F5	VS3	CC-422	Overglaze
CC-423	4R2.5	F9	VS5	CC-423	Medium
CC-524	5M1	F6	VS3	CC-524	Overglaze
CC-525	5M2	F7	VS5	CC-525	Medium
CC-526	5M3	F7	VS5	CC-526	Medium

* CC-000 is keyed to DENTSPLY illuminé bleached shade i7.

**i7 is not a Vita Shade.

± Enamel Shade Stains are optional. Overglaze is required.

¹Vita® and 3D-Master® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. ²Data on file.

CERAMCO IC FIRING CHART

Note: MMT Furnaces, Lower high firing temps by 10°C or 18°F. See firing chart on page 25 if using an Ivoclar™ Programat EP 3000 Furnace.

OPTIONAL: For Ceramco iC All-Ceramic restorations, fill inside of the restoration with heat resistant refractory material (SuperPeg™ II Refractory Material available from DENTSPLY Prosthetics), then insert DENTSPLY Prosthetics All-Ceramic firing pins, and place the restoration on a honeycomb sagger tray to fire. DENTSPLY Prosthetics pins are recommended since they are easy to bend and do not act as a heat sink. Other firing pins will not allow the proper heat transfer and will create a stress point. While initial results with some firing pins may appear acceptable, internal stress can compromise long term success.*

CERAMCO IC FIRING CHART (TO BE USED WITH ALL CERAMCO AND NEY FURNACES)

NP (CoCr) users see chart on page -24

FIRING PROGRAM IN °C

	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Paste Opaque	5 min	5 min	500	55	960	500	960	0	2 min	0
Dentin 1st	5 min	3 min	500	100	840	500	840	30 sec	30 sec	0
Dentin 2nd	3 min	3 min	500	100	830	500	830	0	30 sec	0
Add-On (Dentin/Tissue)	5 min	3 min	500	100	815	500	815	0	0	0
Porcelain Fused-to-Metal (PFM) Margin	5 min	5 min	600	90	880	600	880	0	30 sec	0
Natural Glaze	3 min	3 min	500	100	820	0	0	0	30 sec	0
Press-to-Metal (PTM) Margin	5 min	5 min	600	90	870	600	870	0	0	0
Shade Stain/Overglaze/ Final Margin (PFM)	3 min	3 min	500	100	805	0	0	0	30 sec	0
Shade Stain/Overglaze/Final Margin (All-Ceramic/PTM)	3 min	3 min	500	35	750	500	750	0	1.5 min	0
Enamel cutback(All-Ceramic/PTM)	5 min	3 min	500	35	780	500	780	30 sec	30 sec	0

FIRING PROGRAM IN °F

	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Paste Opaque	5 min	5 min	932	99	1760	932	1760	0	2 min	0
Dentin 1st	5 min	3 min	932	180	1544	932	1544	30 sec	30 sec	0
Dentin 2nd	3 min	3 min	932	180	1526	932	1526	0	30 sec	0
Add-On (Dentin/Tissue)	5 min	3 min	932	180	1499	932	1499	0	0	0
Porcelain Fused-to-Metal (PFM) Margin	5 min	5 min	1112	162	1616	1112	1616	0	30 sec	0
Natural Glaze	3 min	3 min	932	180	1508	0	0	0	30 sec	0
Press-to-Metal (PTM) Margin	5 min	5 min	1112	162	1598	1112	1598	0	0	0
Shade Stain/Overglaze/ Final Margin (PFM)	3 min	3 min	932	180	1481	0	0	0	30 sec	0
Shade Stain/Overglaze/Final Margin (All-Ceramic/PTM)	3 min	3 min	932	63	1382	932	1382	0	1.5 min	0
Enamel cutback(All-Ceramic/PTM)	5 min	3 min	932	63	1436	932	1436	30 sec	30 sec	0

*SuperPeg™II is a trademark of Harvest Dental

CERAMCO IC FIRING CHART (TO BE USED WITH ALL OTHER FURNACE BRANDS. SEE EXCEPTION OF IVOCLAR EP 3000 FURNACE ON FOLLOWING PAGE)

NP (CoCr) users see chart on page -24

FIRING PROGRAM IN °C

	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Paste Opaque	5 min	5 min	500	55	960	500	960	0	2 min	0
Dentin 1st	5 min	3 min	400	60	835	500	835	30 sec	15 sec	0
Dentin 2nd	3 min	3 min	400	60	825	500	825	0	15 sec	0
Add-On (Dentin/Tissue)	5 min	3 min	400	60	810	500	810	0	0	0
Porcelain Fused-to-Metal (PFM) Margin	5 min	5 min	600	90	880	600	880	0	30 sec	0
Natural Glaze	3 min	3 min	500	100	820	0	0	0	30 sec	0
Press-to-Metal (PTM) Margin	5 min	5 min	600	90	870	600	870	0	0	0
Shade Stain/Overglaze/ Final Margin (PFM)	3 min	3 min	500	100	805	0	0	0	30 sec	0
Shade Stain/Overglaze/Final Margin (All-Ceramic/PTM)	3 min	3 min	500	35	750	500	750	0	1.5 min	0
Enamel cutback(All-Ceramic/PTM)	5 min	3 min	500	35	780	500	780	30 sec	30 sec	0

FIRING PROGRAM IN °F

	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Paste Opaque	5 min	5 min	932	99	1760	932	1760	0	2 min	0
Dentin 1st	5 min	3 min	752	108	1535	932	1535	30 sec	15 sec	0
Dentin 2nd	3 min	3 min	752	108	1517	932	1517	0	15 sec	0
Add-On (Dentin/Tissue)	5 min	3 min	752	108	1490	932	1490	0	0	0
Porcelain Fused-to-Metal (PFM) Margin	5 min	5 min	1112	162	1616	1112	1616	0	30 sec	0
Natural Glaze	3 min	3 min	932	180	1508	0	0	0	30 sec	0
Press-to-Metal (PTM) Margin	5 min	5 min	1112	162	1598	1112	1598	0	0	0
Shade Stain/Overglaze/ Final Margin (PFM)	3 min	3 min	932	180	1481	0	0	0	30 sec	0
Shade Stain/Overglaze/Final Margin (All-Ceramic/PTM)	3 min	3 min	932	63	1382	932	1382	0	1.5 min	0
Enamel cutback(All-Ceramic/PTM)	5 min	3 min	932	63	1436	932	1436	30 sec	30 sec	0

Notes: Use the original and modified firing charts for QEX, CeramFire and MMT furnaces.

Use the listed charts for dentin (1st and 2nd) and Add-on for all other furnaces (no hold time for EP 3000 furnace).

Use PTM/AC shade Stain/Overglaze/Final Margin and Enamel Cutback for all furnaces - comments regarding increasing firing temp. and/or time can be found in DFU.

CERAMCO IC FIRING CHART FOR IVOCLAR™* PROGRAMAT EP 3000 FURNACE

FIRING PROGRAM IN °C

	B	S	t[∇]	T	H	V1	V2	L	tL
	Stand-by	Closing Time	Temperature Increase	Holding Temp	Holding Time	Vacuum On	Vacuum Off	Long-Term Cooling	Cool down gradient
Paste Opaque	403	10 min	55	960	2 min	500	959	0	0
Dentin 1st	403	8 min	90	835	0	450	834	0	0
Dentin 2nd	403	6 min	90	825	0	450	824	0	0
Add-On (Dentin/Tissue)	403	8 min	90	810	0	450	809	0	0
Porcelain Fused-to-Metal (PFM) Margin	403	10 min	90	880	30 sec	600	879	0	0
Natural Glaze	403	6 min	90	815	30 sec	0	0	0	0
Press-to-Metal (PTM) Margin	403	10 min	90	870	30 sec	600	869	0	0
Shade Stain/Overglaze/ Final Margin (PFM)	403	6 min	90	800	30 sec	0	0	0	0
Shade Stain/Overglaze/Final Margin (All-Ceramic/PTM)	403	6 min	35	750	1.5 min	450	749	0	0
Enamel cutback (All-Ceramic/PTM)	403	8 min	35	775	1 min	450	774	0	0

FIRING PROGRAM IN °F

	B	S	t[∇]	T	H	V1	V2	L	tL
	Stand-by	Closing Time	Temperature Increase	Holding Temp	Holding Time	Vacuum On	Vacuum Off	Long-Term Cooling	Cool down gradient
Paste Opaque	757	10 min	99	1760	2 min	932	1758	0	0
Dentin 1st	757	8 min	162	1535	0	842	1533	0	0
Dentin 2nd	757	6 min	162	1517	0	842	1515	0	0
Add-On (Dentin/Tissue)	757	8 min	162	1490	0	842	1488	0	0
Porcelain Fused-to-Metal (PFM) Margin	757	10 min	162	1616	30 sec	1112	1614	0	0
Natural Glaze	757	6 min	162	1499	30 sec	0	0	0	0
Press-to-Metal (PTM) Margin	757	10 min	162	1598	30 sec	1112	1596	0	0
Shade Stain/Overglaze/ Final Margin (PFM)	757	6 min	162	1472	30 sec	0	0	0	0
Shade Stain/Overglaze/Final Margin (All-Ceramic/PTM)	757	6 min	63	1382	1.5 min	842	1380	0	0
Enamel cutback (All-Ceramic/PTM)	757	8 min	63	1427	1 min	842	1425	0	0

* Ivoclar™ is not a trademark of DENTSPLY International Inc.

FIRING RECOMMENDATIONS FOR CERAMCO iC - NP (CoCr) AND LASER-SINTERED CoCr**

FIRING PROGRAM IN °C

	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Paste Opaque	5 min	5 min	500	55	960	500	960	0	2 min	0
Dentin 1st	5 min	3 min	500	55	830	500	830	30 sec	30 sec	6 min
Dentin 2nd	3 min	3 min	500	55	820	500	820	0	30 sec	6 min
Add-On (Dentin/Tissue)	5 min	3 min	500	55	805	500	805	0	0	6 min
Porcelain Fused-to-Metal (PFM) Margin	5 min	5 min	600	90	880	600	880	0	30 sec	0
Natural Glaze	3 min	3 min	500	100	810	0	0	0	30 sec	6 min
Press-to-Metal (PTM) Margin	5 min	5 min	600	90	870	600	870	0	0	0
Shade Stain/Overglaze/Final Margin (PFM)	3 min	3 min	500	100	805	0	0	0	30 sec	6 min
Shade Stain/Overglaze/Final Margin (PTM)	3 min	3 min	500	35	750	500	750	0	1.5 min	0
Enamel cutback(PTM)	5 min	3 min	500	35	780	500	780	30 sec	30 sec	0

FIRING PROGRAM IN °F

	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Paste Opaque	5 min	5 min	932	99	1760	932	1760	0	2 min	0
Dentin 1st	5 min	3 min	932	99	1526	932	1526	30 sec	30 sec	6 min
Dentin 2nd	3 min	3 min	932	99	1508	932	1508	0	30 sec	6 min
Add-On (Dentin/Tissue)	5 min	3 min	932	99	1481	932	1481	0	0	6 min
Porcelain Fused-to-Metal (PFM) Margin	5 min	5 min	1112	162	1616	1112	1616	0	30 sec	0
Natural Glaze	3 min	3 min	932	180	1490	0	0	0	30 sec	6 min
Press-to-Metal (PTM) Margin	5 min	5 min	1112	162	1598	1112	1598	0	0	0
Shade Stain/Overglaze/Final Margin (PFM)	3 min	3 min	932	180	1481	0	0	0	30 sec	6 min
Shade Stain/Overglaze/Final Margin (PTM)	3 min	3 min	932	63	1382	932	1382	0	1.5 min	0
Enamel cutback(PTM)	5 min	3 min	932	63	1436	932	1436	30 sec	30 sec	0

** Recommendation: For Compartis CoCr, use NE-Bonder at 980°C

REFRACTORY LAMINATE VENEERS °C

	B Stand-by	S Closing Time	t ^Δ Temperature Increase	T Holding Temp	H Holding Time	V1 Vacuum On	V2 Vacuum Off	L Long-Term Cooling	tL Cool down gradient
Dentin 1st	403	11 min	90	835	2 min	450	834	0	0
Dentin 2nd	403	11 min	90	825	2 min	450	824	0	0

REFRACTORY LAMINATE VENEERS °F

	B Stand-by	S Closing Time	t ^Δ Temperature Increase	T Holding Temp	H Holding Time	V1 Vacuum On	V2 Vacuum Off	L Long-Term Cooling	tL Cool down gradient
Dentin 1st	757	11 min	162	1535	2 min	842	1533	0	0
Dentin 2nd	757	11 min	162	1517	2 min	842	1515	0	0

CERAMCO IC PRESSING CHART

	B	t ^Δ	T	H	E
°C	700	60	890	20 min	100
°F	1292	108	1634	20 min	100

REFRACTORY LAMINATE VENEERS °C

	Pre-dry (Muffle Open)	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Dentin 1st	3 min	5 min	3 min	500	100	840	500	840	0	2min	0
Dentin 2nd	3 min	5 min	3 min	500	100	830	500	830	0	2min	0
Natural Glaze	0 min	3 min	3 min	500	100	820	0	0	0	30 sec	0
Overglaze	0 min	3 min	3 min	500	100	805	0	0	0	2min	0

REFRACTORY LAMINATE VENEERS °F

	Pre-dry (Muffle Open)	Dry (Closing time/min)	Preheat (Closed/min)	Low Temp	Heat Rate	High Temp	Vac Start Temp	Vac Stop Temp	Hold Time Vac	Hold Time Air	Cool Time
Dentin 1st	3 min	5 min	3 min	932	180	1544	932	1544	0	2min	0
Dentin 2nd	3 min	5 min	3 min	932	180	1526	932	1526	0	2min	0
Natural Glaze	0 min	3 min	3 min	932	180	1508	0	0	0	30 sec	0
Overglaze	0 min	3 min	3 min	932	180	1481	0	0	0	2min	0

CERAMCO IC PRESSING REFERENCE CHART

Note: Do NOT preheat ingots

PRESSING REFERENCE CHART IN °C

	Low Temp	*High Temp	Heat Rate	High Temp Hold	Pressing Time 100g Ring Bench Cool	Pressing Time 200g Ring Bench Cool	Repress Time	Pressure	Vacuum
Finesse Press	700	890	60	20 min	20	20	Na	4.25 bars	Full
Multimat Touch & Press	700	890	60	20 min	20	20	Na	2.7 bars	Full
IPS Empress EP 500	700	890	60	20 min	Na	Na	3	5.0 bars	Full
IPS Empress EP 600	500	910	60	20 min	Na	Na	Stop Speed: 300µm/min		Full
QEX CeramPress	700	890	60	20 min	20	20	Na	4.25 bars	Full
Whip Mix Pro-Press 100	700	890	60	20 min	20	20	Na	4.25 bars	Full
CeramFire P	700	890	60	20 min	20	20	Na	2.0 bars	Full

PRESSING REFERENCE CHART °F

	Low Temp	*High Temp	Heat Rate	High Temp Hold	Pressing Time 100g Ring Bench Cool	Pressing Time 200g Ring Bench Cool	Repress Time	Pressure	Vacuum
Finesse Press	1292	1634	108	20 min	20	20	Na	4.25 bars	Full
Multimat Touch & Press	1292	1634	108	20 min	20	20	Na	2.7 bars	Full
IPS Empress EP 500	1292	1634	108	20 min	Na	Na	3	5.0 bars	Full
IPS Empress EP 600	932	1669	108	20 min	Na	Na	Stop Speed: 300µm/min		Full
QEX CeramPress	1292	1634	108	20 min	20	20	Na	4.25 bars	Full
Whip Mix Pro-Press 100	1292	1634	108	20 min	20	20	Na	4.25 bars	Full
CeramFire P	1292	1634	108	20 min	20	20	Na	2.0 bars	Full

* Due to variance in furnaces, high temperatures may need to be adjusted slightly.

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* Ivoclar™ no es una marca registrada de DENTSPLY International, Inc.