GC FUJIVEST SUPER[®]

OCTOBER 2012

R3514

EN Phosphate-bonded investment

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
The logo for GC, consisting of the letters 'GC' in a bold, italicized, sans-serif font, with a horizontal line through the middle of the letters.

DIRECTIONS FOR USE

GC FUJIVEST SUPER® is a carbon-free phosphate bonded investment for precision castings of precious, semi-precious and Pd-base alloys for crown and bridgework, specially developed for quick-heating procedures.

GC FUJIVEST SUPER® featuring:

- ✓ Possibility for use in a quick-heating method or by a progressive heating-up schedule with no difference in expansion.
- ✓ Setting expansion is changed by liquid dilution. Dilution range from 0% (= pure water) to 100% (= pure liquid).
- ✓ Working time (pouring time): more than 4 minutes.
- ✓ Easy mixing because of optimal wettability.
- ✓ Extremely smooth casting surfaces.
- ✓ Very easy digging out.
- ✓ Perfect fitting castings.
- ✓ Possible to use for castable ceramics.
- ✓ Implantology: refer to special Instructions for use "GC Fujivest Super – Implant Instructions for use".

PHYSICAL DATA (typical data)	100% liquid	100% High Exp. Liquid
Setting expansion	1,18 %	1,18 %
Thermal expansion	1,01 %	1,25 %
Total expansion	2,19 %	2,43 %
Compressive strength	5 MPa	8 MPa
Flow	14 cm	13 cm

ISO/DIS 9694.2 (1995)

DIRECTIONS FOR USE

1. Storage	Store powder and liquid at normal room temperature ($\pm 23^{\circ}\text{C}$).	Storage and working temperature of investment powder and liquid are very important in determining the setting time and expansion, and hence the fit and the surface roughness of the castings. Liquid may freeze when exposed to a temperature below 0°C . Once frozen, the liquid cannot be used anymore.
2. Working temperature	Use at $\pm 23^{\circ}\text{C}$ room temperature (19°C minimum).	Higher working temperatures slightly reduce working time but improve the surface smoothness. Lower temperatures prolong setting time and could result in surface roughness.
3. Preparations before investing	Use Multi Sep as a wax separator, which leaves no oily film on the die surface.	
	Use a <u>1mm thick dry casting liner</u> , Casting Ring Liner. Seal the borders of the liner with a thin layer of Vase-line.	Open or wet liners can change the liquid content of the mixture and thus affects the expansion. Be sure that the total inner surface of the metal ring is covered with casting liner. No direct contact of investment with the metal ring.
	If a wetting agent (e.g. GC Surcast) is used, be sure to totally dry the surface prior to investing.	GC Fujivest Super® is a very fluid investment, which can be used without any wetting agents.
4. Expansion	Basic ideas for liquid concentration. The concentration of the expansion liquid depends on various parameters: <ul style="list-style-type: none"> - storage temperature of the investment (powder and liquid) - mixing time - setting under pressure - the amount of mixture 	Pure liquid gives the maximum setting expansion, meaning the largest casting. Pure water gives the minimum setting expansion, meaning the smallest casting. High expansion liquid also available.

Typical properties

	100 % Water	100 % liquid	100 % H.E. liquid
Setting expansion	0,32	1,18	1,18 %
Thermal expansion	0,58	1,01	1,25 %
Total expansion	0,90	2,19	2,43 %

Measured at 23°C

Liquid dilution chart (based on standard liquid)

Pattern	Alloy	Ring size	High precious > 70 % Au		Semi precious < 55 % Au		Pd-base		High precious ceramic alloy		Semi precious ceramic alloy		Pd-base ceramic alloy				
			ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water	ml liq.	ml dist. water			
			x1	x3	x6	x9	x1	x3	x6	x9	x1	x3	x6	x9	x1	x3	x6
1. Inlays/onlays Partial Crowns	- Wax pattern	x1	6,6	6,6	6	7,2	9,2	4									
		x3	17	16	15	18	23	10									
		x6	33	33	30	36	46	20									
	- GC Pattern Resin	x1	8	5,2	7,2	6	10,2	3									
		x3	20	13	18	15	26	7									
		x6	40	26	36	30	53	13									
2. Crown & bridges primary crown	- Wax pattern	x1	7,2	6	6,6	6,6	10	3,2	8	5,2	8,6	4,6	10	3,2			
		x3	18	15	16,5	16,5	25	8	20	13	21,5	11,5	25	8			
		x6	36	30	33	33	49	17	39	27	43	23	50	16			
	- GC Pattern Resin	x1	51	41,4	46,2	46,2	69,4	23	54,4	38	60	32,4	69,4	23			
		x3	8,6	4,6	8	5,2	10,2	3									
		x6	21,5	11,5	20	13	26	7									
3. Secondary parts	- GC Pattern Resin	x1	43	23	40	26	53	13									
		x3	60	32,4	55,4	37	74	18,4									
		x6	7,2	6	7,2	6	10	3,2	7,2	6	9,2	4	9,2	4			
		x9	19	14	18	15	25	8	18	15	23	10	23	10			

Note: All information included in these instructions for use is based on extensive testing and a whole series of casting trials. However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.) different end-results may be obtained.

5. Powder/liquid ratio	Ringsize	Powder	Liquid	Standard P/L ratio = 100gr./22 ml Exact powder/liquid measurement is necessary to obtain stable results. Use electric balance for powder measurement and plastic syringes for liquid measurement.												
	x 1	60 gr	13,2 ml													
	x 3	150 gr	33,0 ml													
	x 6	300 gr	66,0 ml													
	x 9	420 gr	92,4 ml													
6. Mixing	Pre-mix powder and liquid by hand with a spatula.			Be sure all powder is wetted out by the liquid to give a uniform mixture.												
	Place the mixture under vacuum for 15 seconds without mixing.			To reach optimal degassing.												
	Mix for 60 seconds under vacuum (420 rpm).			Always check vacuum level. Insufficient vacuum leads to differences in fit and air-bubbles in the castings.												
7. Working time	At least 4 minutes pouring time at room temperature ($\pm 23^{\circ}\text{C}$).			Working time/pouring time depends on the temperature of powder & liquid and on room temperature. Higher temperatures shorten the working time.												
8. Ringsizes	To use with metal ring. Ringless system only with step heating method. Possible to use with: ringsize x 1 ringsize x 3 ringsize x 6 ringsize x 9			With ringsize x 1 care should be taken. Less exothermic reaction during setting, smaller castings. Temperatures below 23°C may result in rough surfaces. Ringsize x 3 gives the best results. 3 x ringsize = 1 layer of Casting Ring Liner 6 x ringsize = 2 layers of Casting Ring Liner 9 x ringsize = 2 layers of Casting Ring Liner												
9. Investing	Invest under low vibration.			GC FUJIVEST SUPER [®] is very fluid, so strong vibration is not necessary and not advisable. From the moment the ring is totally filled, stop vibration immediately and do not touch the investment until set.												
10. Setting time	Leave to set for 20 min from start of mixing.			Longer setting is possible. Best results are obtained by putting immediately into a preheated furnace after 20 min.												
11. Heating up	Scrape the top surface of the investment ring with a sharp knife.			Diamond disks for trimmers are damaged by the quartz and cristobalite particles of the investment. Due to aggressive burning out, do not open the furnace during heating-up. for Au-alloys for ceramic alloys												
	Put into a pre-heated furnace at end-temperature. ⇒ 750°C ⇒ $800-850^{\circ}\text{C}$ Step by step heating-up schedule is possible.															
	<table border="1"> <thead> <tr> <th>Step heating</th> <th>Heat rate</th> </tr> </thead> <tbody> <tr> <td>1. Room temp. (23°C) → 260°C</td> <td>$2^{\circ}\text{C}-4^{\circ}\text{F}/\text{Min.}$</td> </tr> <tr> <td>2. Holding time at 260°C</td> <td>60-90 min.</td> </tr> <tr> <td>3. 260°C → 580°C</td> <td>$3^{\circ}\text{C}-5^{\circ}\text{F}/\text{Min.}$</td> </tr> <tr> <td>4. Holding time at 580°C</td> <td>30-60 min.</td> </tr> <tr> <td>5. 580°C → 750°C Au-alloy 850°C Ceramic alloy</td> <td>$5^{\circ}\text{C}-9^{\circ}\text{F}/\text{Min.}$</td> </tr> <tr> <td>6. Holding time at end temp.</td> <td>30-60 min.</td> </tr> </tbody> </table> <p>In case of vacuum pressure casting, raise the end temperature by 50°C.</p>				Step heating	Heat rate	1. Room temp. (23°C) → 260°C	$2^{\circ}\text{C}-4^{\circ}\text{F}/\text{Min.}$	2. Holding time at 260°C	60-90 min.	3. 260°C → 580°C	$3^{\circ}\text{C}-5^{\circ}\text{F}/\text{Min.}$	4. Holding time at 580°C	30-60 min.	5. 580°C → 750°C Au-alloy 850°C Ceramic alloy	$5^{\circ}\text{C}-9^{\circ}\text{F}/\text{Min.}$
Step heating	Heat rate															
1. Room temp. (23°C) → 260°C	$2^{\circ}\text{C}-4^{\circ}\text{F}/\text{Min.}$															
2. Holding time at 260°C	60-90 min.															
3. 260°C → 580°C	$3^{\circ}\text{C}-5^{\circ}\text{F}/\text{Min.}$															
4. Holding time at 580°C	30-60 min.															
5. 580°C → 750°C Au-alloy 850°C Ceramic alloy	$5^{\circ}\text{C}-9^{\circ}\text{F}/\text{Min.}$															
6. Holding time at end temp.	30-60 min.															
12. Heating time	Ringsize: x 1 x 3 x 6 x 9			By using quick heating method 40 minutes at end temperature 50 minutes at end temperature 60 minutes at end temperature 90 minutes at end temperature When several investment rings are put into the furnace at the same time, prolong the heating period.												

NOTES

- Clean bowl, spatula, etc. carefully after use. The chemical composition of the residual GC Fujivest Super[®] will delay the setting time of gypsum products.
- Store powder & liquid at room temperature ($\pm 23^{\circ}\text{C}$).**

SAFETY RECOMMENDATIONS AND HAZARD WARNINGS

- Investment materials contain quartz.
Do not inhale dust!
Risk of lung damage (Silicosis, lung cancer).
Advise: wear protection mask Type FFP 2 – EN 149:2001.
Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl.
Rinse the empty investment material bag with water before disposal.
- Remove dust from your working place only when it is wet.
- To avoid the formation of dust when removing the investment material from the casting ring, place the cooled casting ring into water for a short time.
- When sandblasting the cast object, always use a fine-dust filter extraction system.

PACKAGES

- Powder:** 6 kg box (60 gr. Pack x 100)
 6 kg box (90 gr. Pack x 67)
 6 kg box (150 gr. Pack x 40)
- Liquid:** 900 ml bottle
 900 ml bottle High Expansion Liquid

Expiry date: 2 years from the manufacturing

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GC NORDIC AB**Danish Branch**
GC AUSTRIA GmbH**Swiss Office**
GC America Inc.
GC Asia Dental Pte. Ltd.