FINOPRESSURE

Drucktopf

47604

Deutsch

Gebrauchsinformation

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User Information

1. Fields of application

The FINOPRESSURE Pressure Pot is a widely used polymerisation unit in dental technology.

2. Advantages

- Easy handling
- Quick set-up
- Excellent quality of processing

3. Safety precautions

The pot must not be opened when under pressure. Keep out of the reach of children.

4. Processing instructions

The FINOPRESSURE pressure pot is used for cold curing polymerization up to 45 °C. Appropriate for the polymerization of individual trays, occlusal rims, relining, repairs and dentures in slightly warmed water at a max. pressure of 2,2 bars. The pressure pot is made of a special aluminium alloy with a Teflon coating inside and is equipped with an air inlet valve, a safety relief valve and a pressure gauge. Furthermore the FINOPRESSURE has a simple lid closure.

4.1 Commencement of operation

Please check the pot for transportation damages immediately after unpacking. All transportation damages have to be claimed promptly.

Lid and pot handles are not mounted when delivered. When attached with screws, the profiled sides are supposed to face downwards, respectively upwards, where as the even sides should be facing each other. The filling hose is connected to the compressed air supply.

4.2 Handling of the FINOPRESSURE Pressure Pot

The FINOPRESSURE pressure pot is filled with approx. 1 litre of water and the object, which is supposed to be polymerized, is placed inside. Then the pot and the air inlet valve are closed. Screw off the dust protection cap from the filling valve. Now a pressure of 2 – 2,2 bars is set on the pot with the valve head of the filling hose. The pot is only allowed to be heated up to 45 °C.

After the polymerization has finished the inlet valve is opened slowly and the overcharge escapes. The lid may only be taken off the pot, when the pressure gauge reads 0 bars overcharge inside the pot.

4.3 Maintenance

The FINOPRESSURE pressure pot is maintenance free. Lime deposits in the pot can be removed from time to time with lime remover. The pot should not be heated up to more than 45 °C because the Teflon coating could be damaged. It is also advisable not to leave water in the pot over night and to dry the inside if the pot is not used for a longer period of time.

5. Storage

The storage temperature of the unit is determined as -10 °C to +40 °C.

6. Technical data

Max. pressure of compressed air supply: 8,0 bar
Max. operating temperature: 45 °C
Operating pressure: 2,0 to 2,2 bar
Max. allowable operating press.: 2,5 bar
Tested pressure: 4,0 bar
Height: 18,0 cm, with valve 25,0 cm
Max. inside diameter: 20,0 cm
Volume: 4 litres

Weight: 3,77 kg

7. Delivery forms

FINOPRESSURE Pressure Pot with filling hose each 47604

8. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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