FINOSCAN Spray 2
55556/55558

Deutsch

Gebrauchsinformation

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User Information

1. Fields of application

FINOSCAN Spray 2 is a micropulverised scan spray based on well-proven patented spraying technology. It is exclusively for indirect use in laboratories. The spot sprayed on the plaster model is covered with a spraying film. FINOSCAN Spray 2 can be used universally for all CAD/CAM systems and supports an accurate optical scanning quality.

2. Safety precautions

Danger

- H222 Extremely flammable aerosol.
- H229 Pressurized container: may burst if heated.

2.1 Prevention:

- P102 Keep out of reach of children.
- P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.
- P211 Do not spray on an open flame or other ignition source.
- P251 Pressurized container: Do not pierce or burn, even after use.

2.2 Storage:
P410+P412

Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.

Waste codes/waste designations according to EWC/AVV 160504

3. Processing instructions

First put the spraying needle to the aerosol can and turn it into the desired position. Thoroughly clean and dry the surface to be sprayed. Spray from a distance of approx. 5-8 cm. Shake well before use!

4. Storage

Propellant is highly inflammable. Container is under pressure. Protect from direct sunlight and temperatures exceeding 50 °C. Do not force open after use, do not burn. Keep out of reach of children.

5. Technical data

Fill-up quantity:

50 ml net (75 ml gross)

200 ml net (270 ml gross)

6. Delivery forms

FINOSCAN Spray 2

50 ml net (75 ml gross) 55556

200 ml net (270 ml gross) 55558

7. Guarantee

Our technical recommendations of application are based on our own experiences and tests and should only be regarded as guidelines. It rests with the skills and experience of the user to verify that the products supplied by us are suitable for the intended procedures. Our products are undergoing a continuous further development. We reserve the right of changes in construction and composition. It is understood that we guarantee the impeccable quality of our products.

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