Orthocryl®

CE

Gebrauchsanweisung	S. 2
Instructions for use	P. 7
Mode d'emploi	P. 12
Modo de empleo	P. 17
Modalità d'uso	P. 22

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Dear Customer

Thank you for purchasing a Dentaureum product. You have chosen a product that will provide you with consistent high quality, reliability and service.

Optimum performance of our dental technology products is based on working precisely according to the instructions issued by our Research and Development Department. We have therefore included a clear, detailed summary of each working stage in these instructions for use.

If you have any further questions, our team of orthodontic technicians will be pleased to offer advice.

Simply call:

Technical Advice,
fixed appliance technique Hotline Tel. no. +49 (0)7231/803-550
removable appliance technique
Hotline Tel. no. +49 (0)7231/803-555

Information, instructions for use and safety data sheets for Dentaureum products are also available on the Internet at www.dentaureum.com.

Warranty

Dentaureum guarantees a faultless Quality of the products manufactured by us. These non-obligatory recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence in the processing on site. Eventually occurring claims for compensation may only be made up to the value of our goods.

Instructions for use Orthocryl®

1. Use

Orthocryl® Orthodontic Acrylic is a cold-cure resin used in the production of orthopedic appliances.

Optimum stability when the spray-on (salt and pepper) technique is applied. Also suitable for the doughing technique.

2. Composition

Monomer = Methyl Methacrylate „MMA“ (liquid)

Polymer = Polymethyl Methacrylate „PMMA“ (powder)

3. Safety Instructions

See safety data sheets that are supplied with every initial order, or request them at any time.

4. Application

Plaster models should be immersed in warm water for 10 – 15 minutes before applying the resin. This prevents air bubbles from emerging from dried plaster models that can cause the resin to lift from the model, or penetrate the resin.

In order to achieve smooth, silky matt, polished basal surfaces without discolouration or plaster inclusions, we recommend soaking the model in water and then coating with an even layer of the special Orthocryl® separating medium, REF 162-800-00.

Applying the Orthocryl® separator

- After soaking, carefully blow the model dry.
- Pour the separating medium directly from the bottle onto the moist model and spread evenly over the surface using a paintbrush.
- Stand the model onto its back (tuberosities) so that any excess fluid may drain away.
- Carefully blow separation fluid residues or droplets away from the screw and wire retention after approx. 30 seconds have elapsed.

Commence immediately with the acrylic build-up using the spray-on (salt and pepper) technique (A) or the doughing technique (B).

A: Spray-on (salt and pepper) technique

Always start with the powder first, then the liquid, alternating directly onto the prepared model with a gentle pendulum movement.

Important: Only use as much liquid as the powder can absorb. The resin should **not be too wet** to prevent any risk of slumping. Before processing, lightly spray the appliance with powder. This should assure an accurate fit.

B: Doughing technique

Mix the resin in a ratio of 2.5 parts powder to 1 part liquid in a silicone mixing bowl. The well-spatulated mix should be allowed to stand for 5-7 minutes in the covered bowl. Time may be shorter at higher room temperatures. The resin dough is now ready for application.

5. Curing/Polymerization Time

Pressure pot: 2.2 bar, 30 lbs p.s.i.

Water temperature: 40 – 46 °C / 104 – 115 °F

Time: 20 minutes **(no less)**

Polymerisation in a pressure vessel should not be interrupted during the first 12 minutes if possible to prevent the formation of microporosity.

6. Shelf Life

Powder/Polymer: 10 years when container is closed.

Liquid/Monomer: 24 months following date of manufacture in the sealed bottle.

- Store in a cool place at not more than 25 °C / 77 °F.
- Protect from contamination, sunlight or high temperatures as they may cause premature polymerization.

7. Hints

Hints for reducing the remaining monomer content

To reduce the remaining monomer content the completed appliance should be immersed in water for three days, if possible, before the patient uses it. The watering has the effect that the remaining monomer is reduced to values similar to those of heat-curing synthetic resins. Should this not be possible, it is recommended that the patient also places the appliance into water, during the time it is not worn.

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Hints for the hygiene

Appliances should be cleaned regularly and stored in water when not in use to prevent problems with hygiene. This procedure prevents crack corrosion or brown discoloration of the wire components.

Note

Strict quality controls guarantee a uniform quality for all our products.

Orthocryl® Orthodontic Acrylics have been scientifically tested. We have information available on request regarding testing, detailed processing instructions with practical tips, and the entire Orthocryl® product range.

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